## PLACAS CALEFACTORAS DE PRECISIÓN PRECISION HOT PLATES

PLACTRONIC  
PLACTRONIC  
PLACTRONIC

6155100  
6156100  
6157100

# Índice

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# Contents

<b>1 Safety</b> .....	<b>11</b>
2.1 Safety Icons .....	11
2.2 Risks that the operator can be subjected to: .....	11
2.3 Qualified users .....	11
<b>2 General information</b> .....	<b>12</b>
<b>3 Technical specifications.</b> .....	<b>12</b>
<b>4 Contents list</b> .....	<b>13</b>
4.1 Standard contents .....	13
<b>5 Introduction</b> .....	<b>13</b>
<b>6 Equipment description</b> .....	<b>13</b>
6.1 Equipment .....	13
6.2 Features .....	13
6.3 Control panel .....	14
6.4 Rear part .....	14
<b>7 Installation</b> .....	<b>14</b>
<b>8 Operation</b> .....	<b>15</b>
8.1 Starting up .....	15
8.2 Setting an offset temperature .....	15
8.3 Temperature switch .....	15
<b>9 Maintenance</b> .....	<b>16</b>
9.1 Cleaning .....	16
9.2 Spare parts .....	16
<b>10 Guarentee</b> .....	<b>16</b>
<b>11 Electric diagram</b> .....	<b>16</b>
<b>"EC" CONFORMITY DECLARATION</b> .....	<b>18</b>

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
# 1 Safety

The Equipment comes complete with safety features.

This manual indicates areas of possible risk.

## 2.1 Safety Icons

Identifies risk situations and safety measures that should be adhered to. The icons make reference within the paragraph, marked as a lined grey.

### Risk danger

Risk danger

Beware, follow instructions as indicated.

### Electrical risk

Risk of electrical shock if accessing zones are shown with this icon.

Beware, follow instructions as indicated.

### Risk of burns in high temperature areas.

The temperature in the zone indicated with this icon can exceed 60°C. Use thermal protective gloves to complete described functions.

Beware, follow instructions as indicated

### Important information

- Important information in obtaining the best results or optimum performance of the equipment.
- Important information to extend the life of the equipment and maintain its optimum performance.

## 2.2 Risks that the operator can be subjected to:

- Possibility of touching surfaces that are over 60°C.
- Electrical risk.

## 2.3 Qualified users

This equipment should only be used by personnel that are qualified users.

This equipment should only be used by personnel that have understood this manual or have previous adequate instruction of the use of steam sterilisers.

## 2 General information

- 2.1 Handle the equipment with care. Unpack and check that all items coincide with the delivery note. If you see any discrepancy or damage of any kind, notify your distributor as soon as possible.
- 2.2 Read this manual before using the equipment.
- 2.3 These instructions are an important part of the equipment and should be made available to all users.
- 2.4 If in any doubt or clarity on the use of this instrument, contact your local distributor or the service department of J.P. Selecta, s.a.
- 2.5 **ATTENTION, EQUIPMENT NOT CLEAN AND DISINFECTED SHALL NOT BE REPAIRED.**
- 2.6 Any unauthorised modification, lack of maintenance or deliberate damage to the machine, transgresses the directive 89/655/CEE where the manufacturer will not be responsible for damages that could occur due to undue care and attention.
- 2.7 Do not use the equipment for fluids that can produce vapours that can form inflammable or explosive mixtures.

## 3 Technical specifications.

Voltage supply 115-230V 50/60Hz according to the characteristics plate indications.

Código Code	Medidas útiles de la placa cm Plate dimensions cm		Medidas exteriores cm Overall dimensions cm			Consumo Consumption W	Peso Weight Kg
	Ancho Length	Fondo Width	Alto Height	Ancho Length	Fondo Width		
6155100	24	24	13	25	42	700	7.5
6156100	20	40	13	22	57	800	8
6157100	25	60	13	27	77	1000	16.5

## 4 Contents list

### 4.1 Standard contents

The standard equipment consists of the following items:

- |                      |         |
|----------------------|---------|
| - Hot plate:         | 6155100 |
|                      | 6156100 |
|                      | 6157100 |
| - Instruction manual | 80064   |

## 5 Introduction

The rectangular hot plates PLACTRONIC are designed to work with electronic temperature regulation. The circuit allow also working set a working time which starts to discount when the plate reaches to the the set temperature. Specially useful to thermic tests.

The plates are made from grinded duralumin o teflon coated which a greader protections against chemical etching.

The temperature on the hot palte can reach200°C but, in spite of this, they are not devices useful water boiling uses.

## 6 Equipment description


### 6.1 Equipment

The precision hotpataes Plactronic have the following features:

- Pt100 temperature regulation probe
- Thermocouple type K temperature limitation probe
- Teflon coated duralumine plat
- Backing plated heating element of 700, 800 or 1000W depending of the model

### 6.2 Features

- Temperature regulation from room temperature +5°C to 200°C
- Stability  $\pm 0,5^{\circ}\text{C}$
- Resolución 1°C
- Working cycles from 1 minute to 999 minutes
- Electronic temperature limiter

### 6.3 Control panel

1. Temperature and time display
2. Temperature set indication led
3. Time set indication led
4. Modification set key
5. Increasing value key
6. Decreasing value key
7. Start/Stop key
8. Overtemperature pilot lamp indicator
9. Main switch

### 6.4 Rear part

10. Fuse carrier
11. Power supply cord

## 7 Installation

### 7.1 PLACEMENT

Place the piece of equipment near of a socket keeping in mind the its amperage on a leveled surface enough resistant to support the weight indicated on the table of the section 3 (Technical features), leaving and space of 10 cm by the rear part and both laterals.

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### 7.2 POWER SUPPLY CONEXION

Make sure that the voltage of given by the socket is the indicated on the feature label attached on the rear part .

Do not use the equipment unless it is correctly earthed.


Be careful with the plate surface because its surface may produce burnings even the hotplate be switched off.

If you replace the plug keep in mind the following:

Blue wire:	Neutral
Brown wire:	Phase
Yellow/green wire:	Ground

## 8 Operation

### 8.1 STARTING UP

- Connect the equipment to the power supply
  - Turn on the main switch (9)
- For some seconds the display shows «888» and then «-A-»  
Before to continue keep in mind that if the led (2) or the led (3) blinks, it means that the hotplate is heating. If is the case, stop the process pressing the key (7) Start/Stop unless that is not possible modify the set temperature or time.
- The display shows the actual temperature of the plate if the led (2) lights or «000» if the led (3) lights.

#### 8.1.1 SETTING A TEMPERATURE

After doing that is explained in the section 8.1:

- With the key (5) select led (2) lit.
- Press the key (4) and, keeping it pressed, with the keys (5) and (6) select the set temperature which will be shown on the display (1).

#### 8.1.2 SETTING A CYCLE TIME

After doing that is explained in the section 8.1:

- With the key (6) select led (3) lit.
- Press the key (4) and, keeping it pressed, with the keys (5) and (6) select the set cycle time which will be shown on the display (1).

### 8.2 SETTING AN OFFSET TEMPERATURE

When a calibration is made it is possible set an offset temperature:

- Press the keys (5) and (6) simultaneously for at least 5 seconds.
- The display (1) will show the offset temperature (entre -9 i +9 °C)
- While the display shows the offset temperature, this value can be modified using the keys (5) and (6).
- If any modification is done for 5 seconds, the micro processor save the last offset value shown on the display and get out by it self of the adjust menu.

### 8.3 TEMPERATURE SWITCH

The Plactronic hotplates have an electronic temperature switch which avoid that, in case of main temperature regulator failure, the plate reach dangerous temperatures. When the temperature switch triggers the red pilot lamp (8) lights up and cuts the heater. If the problem persists call the technical assistance.

**Do not touch the teflon coated plate or the nearby steel parts. There is risk of burnings**

## 9 Maintenance

Before of remove the cover to handle inside, unplug the the equipment from the socket.

The handle of the electronic circuit boards by unskilled personel may cause damages. Make sure to give the equipment to an authorised technical assistance.

### 9.1 CLEANING

To clean the different parts of the hotplatates we recommend the following products:

- Stainless steel: Alcohol.
- Control panel and teclon coated zone: Alcohol with a cotton clot

### 9.2 SPARE PARTS

The most frequent spares are the following:

Spare description	6155100	6156100	6157100
Main swith	20100	20100	20100
Pilot lamp	22004	22004	22004
Electronic temperature switch	29156	29156	29156
Control circuit board	29327	29327	29327
Heating element	39143	29140	39900
Temperature switch termomouple	43038	43038	43038
Temperature control Pt100 probe	43073	43073	43073
Mains filter	45017	45015	45015


## 10 Guarentee

This product is guaranteed for one year. The guarantee does not cover damage caused by incorrect use or causes beyond the control of J.P. SELECTA,S.A.

Any manipulation of the apparatus by personnel not authorized by J.P. SELECTA,S.A. automatically cancels the guarantee.

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