READER M-2000

4120013

4120014

**LECTOR DE TIRAS MICROTITER
MICROTITER READER**

ÍNDICE

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CONTENTS

- 0. General information
- 1. Packing list
- 2. Basic concepts
- 3. Sample preparation
- 4. Reagents
- 5. Pipettes handling
- 6. Analysis techniques
- 7. Photometer description
- 8. Instalation
- 9. Initial configuration
- 10. Operation
- 11. Warranty

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0. GENERAL INFORMATION

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise your distributor immediately.
- 2) Do not install or use the equipment without reading this handbook beforehand.
- 3) This handbook must always be accessible to all users.
- 4) If you have any doubts or enquiries, please contact your supplier or J.P. Selecta's technical service.
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY APPARATUS TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) If any modification, parts missing or fault due to poor maintenance of any part of the equipment by the user transgresses the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.

1. PACKING LIST

Reference / Part list

Reader
Mains supply cord
Strip holder
Instruction manual
80137

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2.-BASIC CONCEPTS

CELL:

Glass or transparent plastic container, with a thickness and depth precisely defined, for the photometric analysis of liquids. For the mensuration it is introduced into the photometer cell holder.

PHOTOMETER:

Electronic instrument that measures the light intensity across the liquid that is in the cell. This attenuation is due to the absorption of the light at certain wave length. The concentration of the solution contained in the cell is measured by means of this light attenuation.

WAVE LENGTH:

The light spreads in form of periodic waves. The wave length give the colour of the light and is measured in nm (nano-meters). The visible light has a wave length among 380 and 780 nm that it includes all the colors of the rainbow being you blue in the shortest and red in the longest waves. Below 380 nm, we are in ultraviolet region and over the 780nm we are in the infra-red region.

SAMPLE:

Small representative quantities, extracted of the organism or of their excretions to carry out biochemical exams, hematológicos or citológicos.

PLASMA:

Watery residual coming from the separation of the cellular components of the not coagulated blood, for example, for centrifugation.

SERUM:

It is the watery residual of the coagulated blood. The separation of the clot of the serum can be made by decantation or centrifugation.

ENZYMES:

They are proteins that catalize specific biochemical reactions. We measure the enzymatic activity under constant and exactly defined maintained reaction conditions. The substrate should have a high concentration, the pH it should be the good one for the reaction and the temperature should be preset and controlled strictly, the enzymes show a much bigger activity to 37 °C than to 30 °C. According to the type of enzymes identified in the fluids we can deduce the type of affected cell.

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ELECTROLYTES:

They are substances hidrosolubles that transport electric loads to the free state. Cationes have positive load while the anions load negative. Examples: Ca +2, K +, Na+, Cl -.

REACTION TIME

It is the necessary time so that the whole substance to measure is transformed chemically in an appraisable substance with the photometer. For an enzyme: necessary time so that the enzyme can made measurable a quantity of sample.

INCUBATION TIME

It is the necessary time so that the chemical reaction starts in the kinetic techniques

ABSORTION TIME

Time along when the chemical reaction is taking place in technical kinetic.

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3. - PREPARATION OF THE SAMPLE

The blood is the corporal fluid more used in the laboratory with analytic purposes.

For the obtaining of blood, veined punción should be made, avoiding the hemólisis, using needles of the biggest possible caliber and carrying out the extraction slowly. The blood is mixed with anticoagulant shaking it slightly. In a veined extraction, the hemólisis can frequently be caused by these causes:

- Forzing the the blood to cross the tube.
- Shaking the tube too vigorously.
- Extracting blood from an hematoma.

If the sample becomes hemolized interferences dramtically in the hematocrits results, eritrocitos, cholesterol, the enzymes, the iron, etc.

3.1. BLOOD CONSERVATION

Once extracted and blended with the anticoagulant it is left at temperature ambient for an 1 hour, and then stored on to refrigerator at 4°C if it is necessary. The blood should not be refrigerated that when it will be used to obtain serum or plasma.

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3.2. OBTAINING SERUM

It is obtained from a sample of veined punción blood, avoiding the hemólisis and deposited in a recipient without anticoagulant. Left it to resting one hour at ambient temperature. For centrifugación the serum is obtained. The serum can be refrigerated at 4°C, or to freeze at -20°C for a long term conservation. The recongelación of the samples is not recommended.

3.3. OBTAINING PLASMA

It is obtained from the plasma centrifugación of samples of blood with anticoagulant. The plasma with EDTA is more commonly used for hormones, the plasma citrate for tests of clotting.

The plasma heparinized is used for almost all the biochemical tests, even ions.

The plasma with EDTA, used for the realization of the one

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4. - REAGENTS

4.1. CONSERVATION

Most of the reagents should be kept in refrigerator (between 2°C and 8°C). It should not be freezed if it is not specified in the instructions or in the box. There is a series of reagents that is not necessary that they are conserved in refrigerator, but rather they are sufficiently stable to ambient Temperature, whenever this it doesn't exceed the 25°C, in which case, it is always better their conservation in refrigerator. The reagents that don't need conservation in refrigerator are specified in the boxes.

A bad conservation reduces the date of expiration of the reagents considerably.

4.2. EXPIRATIONS

The reagents are stable at least so far of printed expiration for the maker. If the conservation has not been the appropriate one, this it can be lowered, giving erroneous results.

4.3. CONTAMINATION

The reagents can be easy contaminated by a badly manipulation. The sources of contamination can be:

- a) Wrong use of the micropipette:
- The tips should be changed when change the a reagent.
 - no the same pipette should be used for the samples and for the reagents.
- b) The reagent containers should remain open the smallest possible time and avoiding splashes of other reagents or sample.

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WARNING !!!

The test chemicals have they own instruction manual and must be followed futhermore that this general instructions

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5. - HANDLING OF PIPETTES

5.1. TYPES DE MICROPIPETTES.

- a) According to the volume: They can being of fixed or adjustable volume. They are preferable those of volume seleccionable for the versatility that provides to cover a bigger fan of dosage of samples and reagents.
- b) According to the technique:
- pipettes of displacement of air. They work with disposable tips.

It is indispensable to change the tip every time that you different liquid pipeteen, be already reagents or samples, to avoid contaminations.

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- Piston Pipettes with positive displacement. They are pipettes that have the end cap fixed and therefore it is necessary to have much precaution.

In order to avoid contaminations it is necessary to wash the pipette at least three times with distilled water and reject between two pipette usages

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6. TYPES OF TECHNIQUES AND CALCULATIONS METHODS

6.1. TECHNIQUES TYPES

The main types of techniques that are used to analyze samples are: END POINT, DIFFERENTIAL, KINETICS and MULTIESTANDARD. The three first can be of calculation by factor or standard, with which we can classify them in these groups:

END POINT	FACTOR
END POINT	STANDARD
DIFERENCIAL	FACTOR
DIFERENCIAL	STANDARD
KINETICS	FACTOR
KINETICS	STANDARD
MULTIESTANDARD	

In the "end point" techniques the concentration of a certain sample is obtained multiplying the absorbance at a wavelength by a factor.

Sometimes we know this factor (Calculation by FACTOR), other times we have a standard sample and what we know is its concentration (Calculation by STANDARD)

The used expressions are the following ones:

$$\text{FACTOR} \quad C = \text{ABS}_m \times F.$$

$$\text{STANDARD} \quad C = \frac{\text{ABS}_m}{\text{ABS}_{std}} \times C_{std}$$

In the DIFFERENTIALS techniques the concentration values yields from multiplying the difference between the absorbance of a sample and the one of a standard (blank of sample) by a proportionality factor.

As in the previous case we can know this factor (Calculation by FACTOR) or the concentration a standard one (Calculation by STANDARD)

The applied expressions are the following ones:

$$\text{FACTOR} \quad C = (\text{ABS}_b - \text{ABS}_m) \times F$$

$$\text{STANDARD} \quad C = \frac{\text{ABS}_b - \text{ABS}_m}{\text{ABS}_{bstd} - \text{ABS}_{std}} \times C_{std}$$

C	sample concentration.
ABS _m	sample absorbance.
ABS _b	blank sample absorbance.
F	proportionality factor.
ABS _{std}	standard absorbance.
ABS _{bstd}	standard blank absorbance.
C _{std}	Standard concentration.

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In the KINETIC techniques the concentration value is calculated multiplying the reaction speed by a proportionality factor. This speed depends on the existing concentration on the substance that we want to value. The speed of reaction is determined measuring the absorbance variation.

Also in this case we can know the factor (Calculation by FACTOR) or the standard concentration (Calculation by STANDARD)

The formulae are:

FACTOR $C = D \times F$

STANDARD $C = \frac{D \times C}{D D}$

C concentration
 D sample reaction speed
 D D standard reaction speed

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$ABSstd_n < ABSm < ABSstd_{n+1}$


In the MULTIESTANDARD techniques a calibration curve is created by means of several standards.

To analyze a sample we measure its absorbance and we interpolated both in the straight line of concentration of standards of superior and immediately inferior absorbance.

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7. DESCRIPTION

Technical specification

- Wavelength range: 400 to 700 nm.
- It can use interferential filters 10 nm band pass.
- Supplied with filter of 450nm.
- Absorbance range: -0,3 to 3,5 O.D.
- Photometric accuracy :>1%
- Photometric precision: +/- 1%
- Photometric stability: >0.004 A/h.
- Light source: LED diode.
- Detector: Solid state
- Sample compartment: strips for 12 microtiter wells.
- Alphanumeric LCD display
- Concentration measurement by:
 - End point (by K factor or standard)
 - Differential (by K factor or standard)
 - Polygonal curve


General description

The READER M2000 analyzer is a fully-programmable, microprocessor controlled, filter photometer.

Measuring is based on a white light source from LED diode. The monochromator uses interferential 450nm filter and a photodetector with a limited spectral response between 400 and 700nm.

The signal captured by the detector is amplified and filtered to eliminate any interference.

Once the signal proportional to the rate of absorption has been obtained, a digital analogue converter introduces the information into the microprocessor that converts the signal to a logarithmic scale according to Beer's law, that determines the light Absorbance of the sample being measured.

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

7.1 KEYBOARD

The keyboard allows the user to input alphanumerical information into the instrument. The following table describes the various key functions.

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7.2. STRIP HOLDER

Unit in which the Strip are inserted to perform the test.

The light beam cross through cell to quantify the light absorption. This beam goes from up to down.

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8.- INSTALLATION

Choosing the analyser location:

- Place on a flat surface.
- Near to a power source.
- Avoid expose the device to direct sunlight.
- Avoid continuous vibration areas.
- Avoid damp, warm or dusty areas.

Check that the voltage of the analyser, as shown on the label on the rear panel, is the same as that of the mains electricity supply.

Connect the power cable to the socket at the back of the analyser. Next, plug the other end into the mains electricity supply. Ensure that the device is correctly earthed.

To replace the fuse, unplug the mains cable from the analyser. Pull-out the fuse holder. Replace with a new fuse of the same rating as indicated on the specification label.

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9. INITIAL SET-UP

For the propely operation of the analyser is necessary setup some of its functions. Set this functions using "INTERNAL PARAMETERS" menu.

Date and time:

It should be set the right time and date to work with the Photometer.

Language:

We should choose the language in the one which are shown the different messages in the screen.

ABS correction:

It is a factor which affect all the ABS readings. Factory default this factor has a value of 1.00. We can set it with values between 0.80 and 1.20 if measurements with suggest it.

Modify only if your sure.

Measure units:

In addition to 8 measuring units that the analyser incorporates. Two extra units must be added. Use the menu INTERNAL PARAMETERS / MEASURE UNITS.

Printer

Select data send by serial port (RS232) To PC or printer

10. OPERATION

On starting the analyser the DIRECT READING screen appears.It is advisable let the equipment under way during 10 minutes before beginning to work.

Using the keyboard moves for the different menus and submenús.

We make the analysis using a **programmed technique** or use direct readings of ABS and calculate the results manually.

In direct reading the reading of ABS is shown in real time (each second are modernized)

To work with a programmed technique this it should be previously introduced. To see menu PROGRAMMING.

The menu of PROGRAMMING is also good to consult or to modify a technique programmed previously.

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

10.1. MAIN MENU

Use the keys   to move along the menu options.

DIRECT READING
PROGRAMMED READING
PROGRAMMING
MEMORY
INTERNAL PARAMETERS

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10.2 DIRECT READING

In the screen of direct readout is to the babsorbancia that is updated once per second.


On this screen press ENTER. The equipment asks for the blank and next it make the absorbance readings.

This screen is useful to verify the stability of the readings and the nonexistence of drift.

Perform the ABS reading

- Sample must be available in a strip
- Insert the blank.
- Insert the sample.

The ABS reading is shown on the display.

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10.3 PROGRAMMED READINGS

The programmed techniques are used to make analysis whose parameters we have been stored previously in the photometer memory. This way after analyzing a sample the equipment makes all the necessary calculations and in the screen it appears the final result of the analysis (the concentration of a determined substance).

The programmable parameters in each technique are:

1. -Number of technique

It is the position of the memory in which the parameters were stored.

2. -Name

The name that we will assign to technique (up to 9 alphanumeric characters and symbols)

3. -Type of technique

Choose between techniques of END POINT, DIFFERENTIAL or MULTIESTANDARD

4. -Way of calculation

By FACTOR a proportionality factor is introduced. By STANDARD the concentration of a standard is introduced and when selecting the technique the equipment asks for absorbance of the standard.

5.- Duplicate

It make two readings and show the average

6. - Units

They are the units that are in the final result of an analysis. 10 units could be used 8 are preset and 2 can be introduced by the user.

7. - Positive control

It is the upper acceptable value in the technique. It does not appear in the screen when we executed a technique but if it appears in the printed report.

8. - Negative control

It is the lower acceptable value in a technique. It does not appear in the screen when we executed a technique but if it appears in the report.

9. -Standard / Factor

It is the factor of proportionality in the cases of calculation by factor or the concentration of the standard one in the cases of calculation by standard

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How to perform an analysis using a **technique programmed** previously

1. - Select PROGRAMMED TECHNIQUES and to press ENTER.

In screen appears the technique programmed in the position 1 or NO ITEM if there is any programmed technique.

2. - With the arrow keys we moved from a technique to another one, pressing once to move us to the previous or later technique or keep pressed the key for a fast advance. In the screen it appears the number, the name and the type of technique.

Once located the technique press ENTER.

Then appears the filter that it has assigned for the technique, and PUT BLANK (if the parameter PRINTER is set to "YES" there is little dealy necessary to send the parameters of the technique to the PC or the printer)

3. - Place now (one of the two options)

a strip with distilled water if is technique of reading in front of water.


Nothing if is technique of reading in front of air.

Press ENTER. The equipment makes the blank reading and next it ask to introduce into the strip holder the blank sample or the standard if it is a differential technique or with standard, or the sample for the rest of techniques

4. - Introduce next the strip and press ENTER until it appears PUT SAMPLE to place then the sample and to press ENTER. The result of the analysis was in screen

5 - . When pressing to enter the equipment again it asks for a new sample.

When repeating steps 5 and 6 are made new analyses. In the screen it also appears the number of made analyses.

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Only for techniques with Standard.
Only for kinetic techniques with duplicate

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From the screen PROGRAMMING a new technique can be entered (up to 59) or modified.

The programmable parameters on techniques have already seen in section PROGRAMMED READING (10.3)

We will see now how to program a technique and we will explain what is each one of the parameters that take part in the execution of the technique.

1. - In the screen PROGRAMMING press ENTER. Nº 1 with the name appears with the program of the technique if it has been programmed or square black in the nine positions if not.

2. - With the arrow keys we move along the techniques, pressing once to move us to the previous or later technique or maintaining pressed the key for a fast advance. In the screen it appears the number and the name of the technique. When we have selected the position that we wished to press ENTER

3. - In screen appears the black name (in techniques previously programmed) or pictures. In order to introduce the new name to press the arrow keys. We will jump to the first character. To press the arrow keys to modify the first character again (if we maintain the key pressed the symbols advance of continuous way) When we see the symbol that we want to memorizar in the first character we will press ENTER. We jumped thus to the second character. We must repeat this operation until programming the 9 characters of the name. If we want to assign a name with less than nine characters pressing ENTER after assigning the wished characters they introduce spaces in target. For example: GOT we introduce the characters from the 1 to the 3 and from we pressed here to enter to complete the 9.

4. - In the second line of the screen it appears a type of technique (End point,) with the arrow keys we choose the type of technique that we want to program. Press ENTER to accept it.

5. - In the second line of the screen it appears the way of calculation (standard or factor). With the arrow keys we choose the calculation way. Press ENTER to accept. (the value of the standard factor or is introduced later)

6. - Then appear the units. With the arrow keys we choose the units with which we want to work. Press ENTER to accept.

7. - It appears positive control. In order to introduce a value press key CLEAR/ZERO he will appear 0,000 in the screen. With the numerical keys introduce the wished value. We pressed ENTER to accept.


8. - It appears negative control In order to introduce a value of press key CLEAR/ZERO he will appear 0,000 in the screen. With the numerical keys to introduce the wished value. We pressed ENTER to accept.

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
9. - It appears standard factor or. In order to introduce a value press key CLEAR/ZERO he will appear 0,000 in the screen. With the numerical keys to introduce the wished value. Press ENTER to accept.

In order to modify a parameter in a programmed technique we must select the technique and press ENTER until it appears in screen the parameter that we wished to modify. Later we will act as if it was being programmed for the first time (to press keys of arrow or CLEAR/ZERO according to is the case). Once modified we will press ENTER to accept. With key ESC we return to the main menu.

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Push ENTER to accept the upper normal limit
or CLEAR to introduce a new value using numeric key.

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ENTER to accept the lower value or
CLEAR to entry a new value.
Use numeric key.

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ENTER to accept the factor (or Standard) or
CLEAR to entry a new value.
Use numeric key.

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10.5. INTERNAL PARAMETERS

To use this menu to fit some parameters of the analyzer.

How to adjust these parameters? and what is each one of them?

In the main menu we choose INTERNAL Parameters and we press ENTER and with the arrow keys we select the parameter that we want to fit. Press ENTER to validate into the adjustment screen.


1. - DATE and HOUR. It takes the date and time from the internal clock of the equipment. It will appear the date and the hour in the second line of the screen. Press the arrow keys to modify each one of entrances (day, month...) and ENTER to accept it. Press ESC to leave without modifying. Time and date must be properly entered to yield a good and precisely reports.

2. - LANGUAGE. Choose the language of the texts that are in screen. With the arrow keys we choose the language and press ENTER to accept.

3. - CORRECCION ABS. Is a correction factor to apply to the reading in order that the absorbance reading match a standard. Press CLEAR/ZERO and introduce the new value with the numeric keypad. Press ENTER to accept. The value must be from 0,8 to 1,2 otherwise in the screen it appears ERROR.

4. - UNITS OF MEASUREMENT. The equipment reserves two positions of memory so that the user can introduce in them two units of measurement different from which incorporates of series. Press CLEAR/ZERO and introduce the different characters with the arrow keys. Press ENTER to accept.

5. - IMPRESSION. Select data send by serial port (RS232) To PC or printer

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
10.6. MEMORY PRINTING


Prints a hard copy of all programmed techniques or all parameters of a selected programmed technique

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10.7 RS-232 SERIAL OUTPUT

By using the RS232 connection, data can be downloaded to an external computer or directly to a serial printer. The readings are automatically sent after each reading.

- Connection cable:

- Transmission parameters:
Baud rate: 2400; Parity: None;
Nº of bits: 8; Bits stop: 1

- Transmitted characters:
After reading, the analyser sends a message, depending of technique type :

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11.- GUARANTEE

This product is guaranteed for one year from date of purchase against faulty workmanship. The guarantee does not cover damage caused by incorrect use or causes beyond the control of J.P. SELECTA, S.A.

Any mistreatment of the apparatus by unauthorised personnel not approved by J.P. SELECTA, S.A. cancels the guarantee automatically.

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