**BOMBA PERISTÁLTICA PR-2003**

**4129050**

**PERISTALTIC PUMP PR-2003**

**4129050**

## 0. INFORMACIÓN GENERAL

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## 0. GENERAL INFORMATION

1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.

2) Do not install or use the equipment without reading this handbook before.

3) This handbook must always be attached to the equipment and it must be available for all users.

4) If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.

**5) IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY APPARATUS TO BE REPAIRED IF IT IS NOT DULY CLEANED.**

6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.

## 1. CONTENTS

Cap.		Pág.
0.	General information	2
1.	Contens	2
2.	Package list	2
3.	Specification	2
4.	Description	3
5.	Installation	3
6.	Control panel	4
7.	Working with pump	4
8.	Control Menu	7
9.	Dosing	8
10.	Flow	9
11.	Calibration	10
12.	Programmed techniques	11
13.	Feeding	14
14.	Internal parameters	14
15.	CE declaration	15
16.	Warranty	15

## 2. PACKING LIST

The standard package have the following components:

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## 4. DESCRIPCIÓN:

*The peristaltic pump PR-2003 can work in two modes:*

- **Dose:** *The pump doses a set volume.  
Unique dose: It doses a single time.  
Repetition: It repeats the dose every time.*
- **Flow:** *The pump doses a volume for time.  
Programmed flow: During a certain time  
Continuous flow: Continually.*

*The pump allows a calibration that consists on measure in a graduate cylinder the dosed volume and introduce it in the pump.*

*The last calibration is automatically stored.*

*The techniques memory (up to 99) stores the parameters of the most usual applications carried out with the pump.*

*A technique is stored with a name of 5 characters and a number automatically assigned.*

*It allows to choose the work precision: 1 or 0.1 ml. This last is only possible to get it with tubes of interior diameter of 2 mm.*

*It has a dosage range from 0.1 to 9999 ml. It can dose sequential up to 9999 dose with a dose interval from 1 to 9 seconds.*

*The speed of the rotor varies between 60 and 400 rpm. With a flow whose variation will depend in each moment of the interior diameter of the tube used.*

## 5. INSTALLATION

***Before connecting the pump to the mains, check that the voltage matches with the indicated in their characteristics plate.***

***Connect the pump in one plug with protective earth.***


***To change the fuses, use only the ones indicated in the characteristic label.***

***Attention: if the internal clock is not in progress the pump doesn't work.***

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## 7. OPERATION

### Choice of the tube:

It is recommended to use special tubes for peristaltic pumps. Silicone tubes or similar ones with a 2 mm wall have a correct behaviour in peristaltic pumps.

Tubes for peristaltic pump available in J.P. SELECTA, s.a. are:

Code:	ØInt	ØExt
8991203	2	4
8991204	4	8
8991205	6	9
8991206	8	12

### Placement of the tube:

- Open the header.
- Place the tube above the rollers.
- Close the header.
- Adjust the tube's gags of fixation.

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**Programming internal parameters. (See section 14)**

Select the «Language» desired (Spanish, English or French).

Select «Date/Time». Adjust the hour and the date. If it is not adjusted, the pump can work in an anomalous way.

**Feeding: (See section 13)**

Initially the hose is full with air. In order that the dosage will be correct, the hose should be air and bubbles free. So, the bomb will be on until there's no air inside.

**Calibration: (See section 11)**

A calibration must be carried out in order to get reliable dosages.

A graduated cylinder is needed to measure the dosed liquid volume. Then, follow the indications given in section 11.

**Operating without programmed techniques:**

It allows the use of the pump in a quick way for dosage or liquid transfer operations.

**Operating with programmed techniques:**

It is important to have a programmed and calibrated pump in the laboratory when operations are usually carried out. That's the reason why «programmed techniques» are used.

Initially, a «programmed technique» should be programmed and calibrated for a certain tube. The tube's interior diameter is recorded, also, as a «programmed Technique» parameter to remember it later on.

If a technique has been calibrated with a tube, it won't work with another one of different interior diameter. Once programmed, it is stored in the pump memory.

The memorized techniques can be recalibrated time to time without the need of reprogramming the whole technique again.

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**Techniques types**

- Dosage
  - A single dose.
  - Dose repetition every certain quantity of seconds.
- Flow:
  - Continuous operation.
  - Operating for a certain time.

*There are two menus to operate with «programmed techniques»:*

- *Selection menu: to choose a technique from the ones stored in the pump.*
- *Programming menu: to insert a new technique, or to modify or recalibrate an existent technique.*

*When we are operating with a programmed technique, we can consult the number and the name of this technique by pressing key 3 (with standing motor).*

**To program a technique (Section 12)**

*Select «4. - PROG. TECN.» in the main menu.*

*In the “programmed techniques” screen, select “programming”. And then select “new technique”. In a sequential way, different parameters will appear (name, technique type,...). Then, input the new values.*

*It is possible to modify one or more parameters of a “programmed technique” without changing the rest. To do this, press the arrow keys from the “programming” screen to visualize the technique we want to modify and select it.*

*Press ENTER until the parameter that we want to modify appears in screen, and then insert the new value and press ENTER till the next parameter to be modified or to the last parameter.*

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
## 8. CONTROL MENU

When switching the pump on, the screen of main menu appears.

Using keys **q** **p** we can move through the different options

1. - DOSIFICATION
2. - FLOW
3. - CALIBRATION
4. - PROGRAMMED TECHNIQUES
5. - FEEDING
6. - INTERNAL PARAMETERS

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
- **Repetition:** The pump doses the selected volume, «N Dos X» times, with an interval of time programmed in «Interval X sg». Every time the key «ENTER» or the footswitch is pressed, the former sequence is repeated.

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
## 10. FLOW

We can choose two ways to generate a flow:

- **Continuous flow:** The pump starts when pressing «ENTER» and stops when pressing «ESC».

- **Programmed flow:** The pump starts when pressing «ENTER» and stops after a programmed elapsed time.

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
# 11. CALIBRATION

Calibration consists of measuring the dosed exact volume by means of a graduated cylinder. Dosing and Flow are separately calibrated. The last dose and flow calibration are automatically stored and recovered.

This calibration is independent to the one specifically carried out for a programmed technique.

**Without a previous calibration, the pump operation is anomalous.**

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## 12. PROGRAMMED TECHNIQUES


We have the possibility to program up to 99 different techniques.

For each technique, we can choose:

- name of 5 characters,
- if it will be a flow or a dose one,
- the tube's interior diameter,
- the precision,
- the speed, etc.

From the «Programming» menu, we modify the parameters of a programmed technique or we program a new one.

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
-Input the operating precision: 0.1 or 1 ml

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


Time between doses (only when repeated) are selected with up&down keys.

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
## 15. FEEDING

Used to eliminate the air inside the tube, filling it with the liquid to be transferred.

Fix the tubes correctly in their origin and end. Select the pump's appropriate direction. Then press the «ENTER» key. When you notice that the whole air has disappeared from the circuit, press any key to stop the pump.

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## 16. INTERNAL PARAMETERS

In this menu the language is selected among «Spanish, English or French» and the internal clock is adjusted.

**If the internal clock is not adjusted, the pump will not work properly.**

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