**Autoclave clase B para esterilización 18/23 Litros**  
**Class «B» 18/23 liters sterilization autoclave**

**AUTESTER ST DRY PV 18L Clase B**  
**AUTESTER ST DRY PV 23L Clase B**

**4001728**  
**4001731**

# Indice

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# Table of contents

<b>1 Safety</b> .....	<b>14</b>
1.1 Safety symbols .....	14
1.1 Risks which the operator is exposed to.....	14
1.2 Personnel training.....	14
1.3 Proper uses of the equipment.....	14
1.4 Modifications.....	14
<b>2 Packing list</b> .....	<b>15</b>
<b>3 Overview</b> .....	<b>15</b>
3.1 Autoclave "Class B" .....	15
3.2 Programs .....	16
3.3 Bowie Dick test.....	16
3.4 Vacuum test.....	16
3.5 Guidelines to get test an efficient sterilization .....	16
<b>4 Installation</b> .....	<b>17</b>
4.1 Suitable location.....	17
4.2 Connection to power supply .....	17
4.3 Filling water reservoir.....	17
4.4 Emptying water .....	17
<b>5 Control panel</b> .....	<b>18</b>
5.1 The LCD Screen.....	18
5.2 Keypad .....	18
<b>6 Operation</b> .....	<b>19</b>
6.1 Select a cycle (Program) .....	19
6.2 Run a cycle (program) .....	19
6.3 Cycle abort (program stop) .....	19
6.4 System set-up .....	20
6.5 Sterilization end.....	20
<b>7 Maintenance operations</b> .....	<b>21</b>
7.1 Cleaning the draining filter.....	21
7.2 Clean water reservoirs.....	20
7.3 Cleaning /changing the door gasket. ....	21
7.4 Door adjusting.....	21
7.5 Fuse replacing.....	21
<b>8 Technical specification</b> .....	<b>22</b>
8.1 Sterilization report. ....	22
8.2 Error table.....	22
<b>9 Warranty</b> .....	<b>23</b>
<b>10 Statement of conformity CE</b> .....	<b>23</b>

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# 1 Safety

Safety use of autoclaves.

In order to be respected through this manual, the risk of danger situations is remarked.

## 1.1 Safety symbols

It identifies the risk and the safety actions to be applied.

The symbols make reference to the paragraph marked with a grey line.

### Hazard of electrical shock

Electric hazard for entering in areas marked with this signal, or for following the manual instruction marked by this icon.

Follow the instructions for the described operation.

### Hazard of burns to be in contact with hot surfaces

Temperature in this area can exceed 60°C. Use heat protected gloves for the described action.

Follow the instructions for the described operation.

### Important information

- Important information to obtain good results or for an optimal operation of the equipment.
- Important information to extend the life of the equipment or to avoid the degradation of some of its components.

## 1.2 Risks which the operator is exposed to

The operator is exposed to the following risky situations:

- Possibility of touching parts to a temperature above to 60°C.
- Possibility of being in contact with hot steam.
- Electrical risk.

## 1.3 Personnel training

This equipment must only be operated by personnel trained and suitable qualified to work in a laboratory and its common hazards.

This equipment must only be operated by personnel who have read and understood this manual or have been qualified on its manipulation and operation.

## 1.4 Proper uses of the equipment

This equipment is for indoor laboratory use, for sterilizing loads as described.

The sterilization assurance needs also proper procedures for cleaning and storing the loads.

The improper use of this unit (not following manual instructions) can endanger the protection prepared by the equipment.

This equipment its not designed and not manufactured to work in an explosive ambient area.

## 1.5 Modifications

The operation modification or the manipulation of the safety systems unauthorized by the manufacturer, could expose the operator to hazards not indicated in this manual.

## 2 Packing list

Components	Fig.	Qty
• Tray	(1)	1
• Power cord	(2)	1
• Main board Fuse	(3)	1
• Drain water pipe	(4)	1
• Try holder	(5)	1
• Sealing ring	(6)	1
• Fuse		2

## 3 Overview

Class "B" autoclaves are widely used for medical applications, odontologist, veterinary and, generally speaking, to sterilize loads that could be in contact with blood, as "tattoo" and piercing centres, beauty centres, hairdresser, etc...

Class "B" autoclaves assures the properly sterilization for different types of standardized loads: solids, wrapped solids, porous, wrapped porous and hollow devices. To get the sterilization Autester ST DRY PV, uses fraction vacuum system to get an efficient air removal from the load and allows a perfect steam penetration.

The AUTESTER ST DRY PV Class B includes the following features:

- Stainless steel sterilization chamber.
- External cover with epoxy coat.
- Steam generator external from the chamber.
- Digital control.
- Vacuum leakage and Bowie-Dick qualification/diagnostic programs.
- LCD graphic display.
- Filtered air incoming.
- Separated clean/used water tanks.
- RS-232 output for printer connection.

To warranty the load sterilization, the proper program must be used as indicated in this instructions manual.

### 3.1 Autoclave "Class B"

Standard loads for autoclave Class B:

- Solids temperature resistant: Metallic tools, tools.
- Solids temperature sensitive: Plastic devices.
- Wrapped solids.
- Porous: textiles
- Wrapped porous.
- Hollow devices, as catheter, tubes and devices.

This autoclave cannot be used for liquid sterilization.

### 3.2 Programs

The autoclave includes 4 sterilization programs and 2 for validation and diagnostics. Also one program for "Prions" processing is included.

### 3.3 Bowie Dick test

Use this test to validate the autoclave for wrapped and porous loads sterilization.

### 3.4 Vacuum test

Use this test to check the vacuum pump and the chamber leakage for a defined time. Vacuum test is recommended to be done everyday before starting the batch sterilization.

Ensure that the vacuum gauge in the upper left corner of the screen reaches at least a value of -0.7 BAR.

## Programs

Step	B Cycle	18min B cycle	Unwrapped cycle	Wrapped 121°C	BDTest	Vacuum test
Sterilizing temperature	134	134	134	121	134	--
Sterilizing pressure	2.10	2.10	2.10	1.10	-	< -0.8
Class	B	B	S	B	B	--
Pre-Heating	0-7 min	0-7 min	0-7 min	0-7 min	0-7 min	--
Vacuum time	10 min	10 min	4 min	10 min	10 min	4 min
Admission time	5 min	5 min	5 min	5 min	5 min	--
Sterilizing time	4 min	18 min	4 min	16 min	3.5 min	Keep P. 10min
Drying time	10 min	10 min	10 min	10 min	1 min	--
Quick drying time	5 min	5 min	5 min	5 min	--	--
Total time	36 (31) min	50 (45) min	30 (25) min	48 (43) min	19.5 min	15 min
Usage (General)	Class B	Prions	Solids loads	Rubber&plastic	OO*	OO*

\*OO: Operation qualification

### 3.5 Guidelines to get test an efficient sterilization

The material to be sterilized needs to be clean, free from any encrusted or residual material, wash it in a good detergent and distilled water, then rinse with distilled water.

Do not overload the trays, left a free space of 1 to 2 cm between goods to allow a free steam circulation and to ease drying.

#### 3.5.1 Sterilization of instruments in bulk

Place the instruments on a tray with water repellent paper leaving a space between each instrument. Try to keep similar type material instruments together on a tray.

#### 3.5.2 Sterilization of instruments in pouches or bags

Place only one instrument per bag.

Place the bags on water repellent paper in a vertical position so that they do not touch each other or the wall of the autoclave. Better use a rack.

#### 3.5.3 Tubing



Tubing need to be clean and drained. Place tubing so that the openings on both ends are open. Make sure that the tubes are not touching the autoclave wall, blocked or twisted and that nothing can block the inside.

#### 3.5.4 Flasks

Never put hermetically sealed flasks in the autoclave.

Place the contents open end down to prevent any residual water being left inside.

## 4 Installation

Once the unit is unpacked, visually check the packing list (See 1.1) and that there are no parts damaged during transportation.

Save the packing for some days for later sending.

### 4.1 Suitable location

Place the unit over a flat level and stable surface, suitable for the unit weight.

Use the levelling instrument to check the flat placement.

**Note:** The front of the machine can not be lower than the rear part. In this case the unit does not operate properly and could display E3 message (Dry not completed).

No longer that 1.5m from a power supply plug.

No tap water and waste will be needed.

Autoclave dimensions are: high 460mm x wide 460mm x depth 600mm.

To work comfortably, let a free space of 50cm around the equipment.

To allow ventilation, let not least than 20cm on the rear part.

Do not block or cover any ventilation holes.

Prepare to refill the water reservoir from the top of the unit.

### 4.2 Connection to power supply

Locate a power supply socket (rated to 16A) near to the unit.

Use the power cord supplied or another of similar specifications.

For safety reasons, the socket must have ground terminal, since the unit must work grounded.

Do not manipulate the power cord with wet hands.

Before plugging on the socket, check that its voltage match the voltage marked on the unit plate.

Plug the power cord on the autoclave rear part.

### 4.3 Filling water reservoir

Use distilled or deionized water with this autoclave.

Open the top cover for the reservoir and fill it.

### 4.4 Emptying water

Draining of water from water tank and used water tank.

Connect the water outlet using the one end of the tube without a connector. Water tank will drain water after connected, see (Fig. 4).

Used water tank will drain water after turning the draining knob open see (Fig 4.4).

## 5 Control panel

The control panel allows function selection and shows the distiller alarms status.

Once the unit is plugged in, put it on, by pressing the main switch. The unit goes to stand by status (ON).

**Note:** If main switch is pushed consecutively with an interval of less than 3 seconds, the unit could start in abnormal mode. So wait for a few seconds to restart the unit.

### 5.1 The LCD Screen

The LCD screen shows information on depending of autoclave status:

Three states are defined:

- Working state: The autoclave is running a cycle.
- Menu state: User is browsing for the menus.
- Ready to work state: Cycle is selected, waiting "START" key.

Message row 1:

- Working state: Shows temperature and pressure data
- Menu state: Shows menu row or page layout name
- Ready to work state: Shows sterilizer cycle name

Message row 2:

- Working state: Shows pressure curve
- Menu state: Shows menu row
- Ready to work state: Shows selected cycle parameter

Message row 3:

- Working state: Shows pressure curve
- Menu state: Shows menu row
- Ready to work state: Shows selected cycle parameter

Message row 4

- Shows key help information

### 5.2 Keypad

Some keys change their function according to the autoclave status.

#### Menu/ C

- Working state: EXIT key.
- Menu state: CANCEL key. Back to previous menu.
- Ready to work state: ENTER key. Validation input.

#### 5.3 OK

- Menu state: Means OK (Next menu).

#### 5.4 START/STOP

- Start key. It will start sterilization cycle.

#### 5.5 Quick

- Ready to work state: Sets "quick dry" on sterilization cycle.
- Menu state: It means down key.

#### 5.6 STEP

- Working state: Jumps to next step.
- Menu state: It means up key.

J.P. SELECTA, s.a.  
  
AUTESTER  
ST DRY PV

SELECT A CYCLE  
UNIVERSAL B CYCLE  
18 MINUTES B CYCLE  
UNWRAPPED CYCLE

SELECT A CYCLE  
RUBB & PLAS.  
CYCLE  
B&D TEST CYCLE  
LEAK TEST CYCLE

SELECT A CYCLE  
CLEAN GENERATOR

NOTICE  
REFILLING

Fig 6.3

## 6 Operation

After pressing the main switch a message shows the autoclave model and automatically jumps to the main menu.

### 6.1 Select a cycle (program)

- Open the door. Place the load trays.
- Press "Menu/C" key, to enter cycle selection screen.
- Press "STEP ▲" or "DOWN ▼" key to select one of the cycles:

**UNIVERSAL B CYCLE:** 134°C 4 min.

- Press OK key for confirmation.
- Set drying time to 5 or 10 min by pressing "QUICK" key.

**18 min B CYCLE:** 134°C 18 min.

- Use for Wrapped load and Prions processing.
- Press OK key for confirmation.
- Set drying time to 5 or 10 min by pressing "QUICK" key.

**Unwrapped load.** 134 °C 4min.

- For not wrapped load thermal resistant.
- Press OK key for confirmation.
- Set drying time to 5 or 10 min by pressing "QUICK" key.

**Low temperature wrapped load.** 121°C 16min.

- For wrapped load not compatible with 134°C temperature.
- Press OK key for confirmation.
- Set drying time to 5 or 10 min by pressing "QUICK" key.

**Bowie Dick.** Test for sterilization validation or OQ stage.

- Press OK key for confirmation.

**Vacuum test.** Test for leakage checking or OQ stage.

- Press OK key for confirmation.

### 6.2 Run a cycle (program)

Once the load is placed into the chamber and program selected.

- Close the door.
- Press "START/STOP". At the moment the electrical lock acts.
- For safety reasons, if the door locking was not right, the programs do not start.

### 6.3 Cycle abort (program stop)

Once the program is initiated can also be stopped.

- Close the door.
- Press "MENU/C" for 3 seconds. This causes program jumps to the following phase.
- Repeat operation up to the end.

#### Note:

If the water level is lower than the minimum level, the message "NOTICE Refilling" appears on the screen and an alarm will sound.

Afterwards, refill the water reservoir. Do not forget to empty the reservoir for used water.

## 6.4 System setup

For system setup, time, date, printer, etc. the autoclave must not be running a program (Menu state)

- Press "MENU/C" for 3 seconds to accede SETUP menu.

Use keys «STEP ▲» o «QUICK ▼» to select among the menu options.

Use keys «Menu/C» o «OK» to reject or validate an option.

- Select «DATE» to enter the date.
- Select «TIME» to enter the time.
- Select «PRINTER» to activate the printer report.
- Select «Time setup» to allow the sterilization time extension from 0 to 3 minutes added.
- Select «Temperature setup» to adjust temperature reading. Value from -2 to +3°C.

### Note:

Sterilization time modification must be properly validated for each processed load.

### Note:

Do not modify the temperature offset if you do not have a reliable means to measure the inner chamber temperature.

Any modification of these two parameters can cause an inefficient or defective load sterilization.


## 6.5 Sterilization end

After sterilization, the autoclave depressurizes the chamber, then you can open the door.

When the chamber does not have dangerous pressure inside (< 0.05 bar) shows this situation on the screen and unblocks the electrical lock.

If the printer is attached and ON, printed report is produced.

Use provided holder to remove the hot trays from the chamber.

**Caution: when unloading material from the autoclave.**

**Hazard of burns to be in contact with hot surfaces.**

Temperature in this area can exceed 100°C.

Use the tray holder to load or unload the tray with instruments in it, in order to avoid burning injuries.

Loosen/Tigh with hand.

Fig 7.1 Draining filter.

Filter components.

Remove 5 screws.


Fig. 7.2 Cleaning reservoir.

Fig. 7.3 Door gasket replacing.

Fig. 7.4 Door adjusting.

Fuses

Fig. 7.5 Fuse replacing.

## 7 Maintenance operations

To get an optimal performance from the autoclave, program the following maintenance operations.

### 7.1 Cleaning the draining filter

Follow these instructions:

- Locate this filter inside the chamber at the bottom (See fig 7.1).
- Clean approx. on 3 month basis.
- Filter can have inlay or dirty, due to the loads processed.
- Unscrew (use only hand no tools) and pull the filter to remove it.
- Remove the blue rubber and clean.

#### Note:

To avoid problems caused by dirtiness:

- Clean the goods before to sterilize it.
- Specially clean free form oil and grease.
- Use only deionised/ distilled water.

### 7.2 Clean water reservoirs.

Follow these instructions:

- Remove the 5 screws from the cover. See Fig 7.2.
- Lift the cover.
- Clean inlays and with a cloth.

### 7.3 Cleaning /changing the door gasket.


Follow these instructions:

- Check, regularly, that it looks clean and no particles added on it.
- Clean with a clean and smooth cloth.
- Replace for a new one if you see any derated part.

### 7.4 Door adjusting.

If you note that there are steam leakages through the door, you can adjust it by doing the following:

- 1 Stop the autoclave and let it cold.
- 2 Insert the spanner into the gap beneath the plastic door cover.
- 3 Turn the door adjusting nut counter clockwise using the supplied spanner as showed in the pictures.

### 7.5 Fuse replacing.

Follow these instructions:

- Unplug the unit from the mains.
- Use a screwdriver to turn 1/4 the fuse holder cover.
- Use fuse 6x30mm 12A.
- Insert the fuse and fix it.

## 8 Technical specification

Model:	4001728	4001731
Mains voltage:	230V 50Hz (See plate) Other voltages could be available	230V 50Hz
Electrical power:	1800W (See plate)	2300W
Maximum consumption:	8.5 A	10 A
Fuses:	T 12A dimensions: 6x30 mm	T 16A
Weight (Net):	46.6 Kg.	48 Kg.
Dimensions: W x L x H	460 x 460 x 600 mm	495 x 550 x 760mm
Chamber volume:	18 Litres	23 Litres
Chamber dimensions	Ø 250mm x 355mm	Ø 250mm x 450mm
Water reservoir volume:	3.5 litres	3.5 litres
Ambient conditions:	Inner use. Altitude up to 2000m. Temperature from 5°C to 40°C.	
Relative maximum damp	80% for above 31°C temperatures, decreasing lineally to 50% for 40°C.	
Overvoltage category:	Category II	
Pollution degree:	2	
Chamber material & Steam generator	Stainless steel.	

### 8.1 Sterilization report

Date:	Date.
C.N.:	Cycle name.
C.S.T.:	Cycle start time.
Vp1 Vp2 Vp31:	Vacuum level.
Pp1 Pp2 Pp31:	Pressure level.
H.S.T.:	Sterilization starting time.
Max.:	Maximum value pressure and temperature recorded.
Min.:	Minimum value pressure and temperature recorded.
Ave.:	Average pressure and temperature.
H.E.T.:	Sterilization starting time.
D.S.T.:	Drying cycle start time.
C.E.T.:	Cycle end time.
Result:	Cycle evaluation.
Operator:	Operator signature.

### 8.2 Error table

Item	Message	Failure	Checking
1	INVALIDATION	Wrong operation	Check operation according to instructions manual
2	CLOSE DOOR	Door is not closed before running	Check if door is properly closed
3	CLEAN REQ.	Clean operation required	Run cleanning cycle. See 6.1
4	REFILLING	No water on reservoir	Refill tank
5	USED WATER	Used water tank fill	Drain used water tank
6	SENSOR ERROR	Temperature > 250°C	Temperature sensor
7	PRESSURE OVER	Over pressure (>2.6bar)	Pressure sensor. Temperature sensor
8	CHA.TEMP.OVER	Chamber Overtemp. > 160°C	Temperature sensor
9	ST. FAULT	Temperature out of range -1 +3°C	Adjust temperature (Ver 6.4)
10	VENTING ERROR	Venting route is blocked	Input air valve
11	DOOR ERROR	Door open during cycle	Door switch
12	STEAM TEMP.OVER	Steam generator Overtemperature	Temperature sensor
13	WARM UP TIME OVER	Pre- heat faillure>20min	Heater wiring or heating elements
14	PRE.TIME OVER	Steam input time >40min	Check input steam pipe

## 9 Warranty

The warranty covers the defective components or faulty assembling of the unit for 1 year since the purchasing date.

Keep purchasing bill as it is the only valid document for the purchasing date.

The warranty does not cover the damages caused for a unit malfunction. The unit must be suitably qualified before being used for valuable samples processing.

Return good are not accepted without previous J.P. SELECTA, s.a. approval.

## 10 Statement of conformity CE

The manufacturer:

J.P. SELECTA, s.a.  
Ctra. NII Km 585.1  
08760 ABRERA (BARCELONA)  
ESPAÑA

Declares that the laboratory equipment:

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Meet the following European Directives:

73/23/CEE	Directive for electrical safety.
89/336/CEE	Directive for electromagnetic compatibility.(EMC)
97/23/CEE	Directive for pressure devices. Paragraph 3.3

Meet the following Standards:

EN 61326:2002  
Electrical equipment for measurement, control and laboratory use. EMC Requirements.

EN 61010-1:2002  
Safety requirements for electrical equipment for measurement, control and laboratory use. Part 1: General requirements.

EN61010-2-041  
Particular requeriments for laboratory autoclaves.

Abreira, April 2007

Ramón Ramón  
Dtor. Técnico.

David Pecanins  
Rble. Calidad

*J.P. SELECTA, s.a.*

---

---