MECHEROS BUNSEN BUNSEN BURNERS

**7000044
7001106
7000134**

**7000937
7001107
7009133**

Indice

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1. DESCRIPTION

BUNSEN burners are widely used like heat source in laboratories. Its application goes from the incineration of samples to the modeled of glass tubes. When correctly used, they offer a great security and ease of use.

There are different types of burners, of similar characteristics, but sometimes they better fit to a specific job.

2. BURNERS TYPES

BUNSEN burners:

| Code | GAS | Ø Tube | Nominal pressure (mbar) | Heat power |
|---------|-------------|----------|-------------------------|------------|
| 7000044 | Butane | Tube Ø11 | 28/37 | 0.69KW |
| 7000937 | Natural gas | Tube Ø11 | 18 | 0.69KW |
| 7000134 | Butane | Tube Ø11 | 28/37 | 0.69KW |
| 7009133 | Natural gas | Tube Ø11 | 18 | 0.69KW |

MEKER burners:

| Code | GAS | Ø Tube | Nominal pressure (mbar) | Heat power |
|---------|-------------|----------|-------------------------|------------|
| 7001106 | Butane | Tube Ø22 | 28/30 | 0.66KW |
| 7001107 | Natural gas | Tube Ø22 | 20 | 0.59KW |

3. INSTALLATION

- DO NOT CONNECT a gas burner without a pressure reduction.
- DO NOT CONNECT a gas burner directly to the gas cylinder.
- Before connecting the burner to a gas installation be sure that it provides the gas and the suitable pressure.
- The connection elements and tubes must fulfill the requirements established by the providing gas company.
- The installation of gas burners must be made by an authorized installer. Consult with the providing company any doubt in the installation.
- Do not modify any piece (orifices, steps, etc) of the gas burners.


In the table of section 2 the pressures of use are indicated.

4. OPERATION

Before using the gas burner, carry out some previous verifications:

- Check that all the burners have their gas regulation key closed before opening the installation passage key.
- Check that the tube of gas entrance is correctly fixed to the gas entrance tube mouth.
- Check that the installation passage keys are open.

4.1 BURNER IGNITION

- Close totally the air intake 2, open slightly the gas passage key 1 and approach, laterally, an ignited match to the mouth of tube 5.
- Regulate the gas key 1, until obtaining a flame with the wished height.
- Gradually, open the air intake 2.
- DO NOT OPEN IT IN A SUDDEN WAY BECAUSE THE BURNER CAN BE EXTINGUISHED.
- In order to obtain a higher temperature, open more the gas key 1 and air intake 2.
- THE BURNER EXTINGUISHES WHEN CLOSING THE GAS KEY 1.
- In models that have no gas regulation key, please use the installation key.

4.2 FLAME ADJUSTMENT

In order to obtain a flame with great calorific power, the suitable air intake must be adjusted to the gas entrance. If the burner burns with the air intake too much closed, the combustion is incomplete and the flame shows an orange color.

When opening the air passage, and when turning piece 2, the combustion is complete and in the flame two clearly separated zones by a light blue cone are appreciated.

Outside the cone, the combustion is complete, and a great calorific power is obtained. Inside the cone, the combustion is not complete and the temperature is inferior.

Once the air intake is adjusted, it can be fixed with the piece 4 counterclockwise.

How to make a practical verification of the flame adjustment:

Hold a porcelain capsule with a clamp upon the flame:

- If the capsule is blackened by the smoke, it indicates that the air intake is insufficient.
- If some small drops of water arise, it indicates that the air intake is correct.

4.3 MAINTENANCE

Approximately once a year, clean up the injector hole and the outlet pipe of the flame. Check the status of the toric joints of the gas passage key.

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