1. Descripción

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
Fig 1. Partes del Mechero.

1. Overview

BUNSEN burners are widely used like heat source in laboratories. Its application goes from the incineration of samples to the modeled of glass tubes. When correctly used, they offer a great security and ease of use.

These BUNSEN burners are with air regulators and needle valve.

2. Gas pressure

Gas inlet pressures for all types of burner:

BUNSEN burners:

Code	GAS	Nominal pressure (mBar)	
1000876	Natural gas	18	to 25
1000875	Butane	47.5	to 57.5

Warning!

Laboratory burners must be kept under constant supervision. Open the main gas tap before use. And after use, immediately close the main gas tap.

3. Installing de burner

- DO NOT CONNECT a gas burner without a pressure reduction.
- DO NOT CONNECT a gas burner directly to the gas cylinder.
- Before connecting the burner to a gas installation be sure that it provides the gas and the suitable pressure.
- The connection elements and tubes must fulfill the requirements established by the providing gas company.
- The installation of gas burners must be made by an authorized installer. Consult with the providing company any doubt in the installation.
- Do not modify any piece (orifices, steps, etc) of the gas burners.

In the table of section 2 the pressures of use are indicated.

4. Operation

Before using the gas burner, connect it to the gas supply.

Close the air inlet (2). Turn the needle valve (3) fully clockwise. Then turn on the gas at the main tap.

If a Safety Burner is being used, press and hold down the red knob. Then ignite the gas. Hold the red knob down for 3-5 seconds after ignition.

The permanent flame burners keeps a little flame ignited, when gas valve is openend the main flame ingines autotomatically.

Adjust the gas flow with the needle valve (Fig 1 part 3).

Adjust the flame by opening or closing the air inlet. (Fig 1 part 2)


Fig 1. Burner parts.

5. Types of flame

In order to obtain a flame with great calorific power, the suitable air intake must be adjusted to the gas entrance.

For a roaring flame, open the air inlet. For a luminous flame, close the air inlet.

6 Average consumption

Pressure variations can affect the average consumption of butane or natural gas.

Butane/propane consumption:	26 L/h	0.0596 Kg/h
Natural gas consumption:	80 L/h	0.0560 Kg/h

7 Maintenance

Approximately once a year, clean up the injector hole and the outlet pipe of the flame. Check the status of the toric joints of the gas passage key.