FOTÓMETRO ANAL. ENOLÓGICO M-3000 4120019
OENOLOGICAL ANALYSER PHOTOMETER

FOTÓMETRO ANALIZADOR 280NM. 4120023
ANALIZER PHOTOMETER 280NM.

COLORÍMETRO DIGITAL CLORMIC 4120009
DIGITAL COLORMETER CLORMIC

COLORÍMETRO DIGITAL SP. PRO-A 520NM 4120029
DIGITAL COLORMETER SP. PRO-A 520NM

COLORÍMETRO DIGITAL SP. 660NM 4120028
DIGITAL COLORMETER SP. 660NM

INDICE

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9-pin Hembra
nº pin:
2
3
5

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Baud rate: 1200;
Parity: None;
Nº of bits: 8;
Bits stop: 1

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0. GENERAL INFORMATION

- 1) *Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise your distributor immediately.*
- 2) *Do not install or use the equipment without reading this handbook before.*
- 3) *This handbook must always be accessible to all users.*
- 4) *If you have any doubts or enquiries, please contact your supplier or J.P. Selecta's technical service.*
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) *If any modification, parts missing or fault due to poor maintenance of any part of the equipment by the user transgresses the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.*

1. PACKING LIST

	Part No.
1. Photometer / colormeter	
2. Mains supply cord	0007001
3. Polystyrene cuvettes (4 units)	5100022
4. Instruction manual	80269

2. DESCRIPTION

General description of the equipment:

The photometer / colormeter is a small benchtop analyser.

Measuring is based on a light source from a transmitter, a monochromator composed of interference filters and a silicon photodetector with a limited spectral response according to technical specifications.

The signal captured by the detector is amplified and filtered to eliminate any interferences, an analogue circuit converts the signal to a logarithmic scale according to the Beer's law that determines the light absorbance of the sample being measured.

Once the signal proportional to the rate of absorption has been obtained, a digital analogue converter introduces the information into the microprocessor.

The software can measure ABS (Absorbance), Transmittance and Concentration values.

General technical specifications:

- *Manual filter wheel.*
- *Absorbance range: -0,3 to 3,5 O.D.*
- *T% reading range: from 0 to 100 T%.*
- *Photometric accuracy: 1%*
- *Photometric precision: $\pm 1\%$*
- *Photometric stability: >0.004 A/h.*
- *Detector: Solid state*
- *Sample compartment: 1 & 10mm cells.*
- *Minimum volume: 1 ml.*
- *Alphanumeric LCD display*

- Concentration measurement by:
Absorbance by K factor.
Absorbance by factor from standard solutions.

Specific technical specifications:

4120009

Wavelength range: 400 to 800 nm by using interference filters with a 30 nm band pass.

Filters: 420, 440, 490, 520, 550, 580, 620 & 680 nm

Light source: Long life tungsten lamp

4120019

Wavelength range: 340 to 750nm, by using interference filters with a 10nm band pass.

Filters: 340, 420, 520, 620 & 750 nm

Light source: Halogen lamp

4120023

Wavelength range: 280 nm

Filters: 280 nm

Light source: UV photodiode

4120029

Wavelength range: 400 to 800 nm by using interference filters with a 30 nm band pass.

Filters: 520 nm

Light source: Long life tungsten lamp


4120028

Wavelength range: 660 nm

Filters: 660 nm


Light source: tungsten lamp

Main components:

- 1- 16 Key keyboard.
- 2- 40 character display.
- 3- Selected filter.*
- 4- Cuvette holder.
- 5- Filter selector control.*
- 6- Main switch.
- 7- Mains connector.
- 8- Fuse and spare fuse.
- 9- RS 232 output.

* Only models 4120009 & 4120019

2.1 KEYBOARD DESCRIPTION

The keyboard allows the user to input alphanumerical information into the instrument. The following table describes the various key functions.

KEY	ACTION
0 al 9	Numerical entry
Escape (Esc)	Exit to main menu.
Enter (Enter)	Validate entry
Clear / cero	Without function
UP / DOWN	Select menu option
Paper / .	Without function

3. INSTALLATION

Place the analyser on a flat surface near a power source.

Do not expose the device into direct sunlight.

Avoid installing the instrument in excessively damp, warm or dusty areas.

Check that the voltage of the analyser, as shown on the label on the rear panel, is the same as that of the mains electricity supply.

Connect the power cable to the socket at the back of the analyser. Then, plug the other end into the mains electricity supply. Ensure that the device is correctly earthed.

To replace the fuse, unplug the mains cable from the analyser. Pull-out the fuse holder. Replace with a new fuse of the same rating as indicated on the specifications label.

4. OPERATION


Sample analysis:

- Samples must be placed in a:
 - Rectangular 10mm cell. OR
 - Polypropylene/glass Ø12-16mm tube.
- Insert the sample in the cell holder.
- Select the desired wavelength filter.
- Select on the main menu the type of result required.

On starting the analyser, ABSORBANCE appears from the main menu.


Options from main menu can be selected through the keyboard.

Prior to starting measuring samples, choose the most adequate filter that corresponds to your samples wavelength.

4.1 MEASURING ABSORBANCE


To measure absorbance, manually select the appropriate filter.

For absorbances higher than 3.5 and lower than -3.5, the display shows --- (Out of range)


4.2 MEASURING TRANSMITTANCE

To measure transmittance, manually select the appropriate filter.


4.3 MEASURING CONCENTRATION BY FACTOR

Calculates the sample concentration by measuring the absorbance and applying a pre-introduced K factor. Preselect the appropriate wavelength filter. To get good precision, use the same type of flask with distilled water to made the blank.


4.4 MEASURING CONCENTRATION BY STANDARDS

Calculates the sample concentration measuring the absorbance and applying a pre-introduced concentration factor. Manually select the appropriate wavelength filter. To get good precision, use the same type of flask with distilled water to made the blank.

4.5 MEASURING ABSORBANCE IN CONTINUOUS READING


Continuous reading of the absorbance value of a sample every 4 times per second. Manually select the appropriate filter.

4.6 ENTERING A CORRECTION FACTOR

This menu allows the user to enter a correction factor that is applied to the absorbance reading. The default value is 1.00. The value needs to be between 0.8 to 1.2.

Entering this factor is specially useful for repetitive samples whereby an adjustment is required by a standard.

4.7 RS-232 SERIAL OUTPUT

By using the RS232 connection, data can be downloaded to an external computer or directly to a serial printer. The readings are automatically sent when the «ENTER» key is pressed after each reading.

- Connection cable:

9-pin Female
pin number:

2

3

5

- Transmission parameters:

Baud rate: 1200;
Parity: None;
N° of bits: 8;
Bits stop: 1

- Transmitted characters:

After reading, the analyser sends a message, depending on the reading type selected:

ABSEX.XXE	10 chars + CR LF	E: space
TRAEX.XXE	10 chars + CR LF	E: space

4.8 PERFORM A RESET

To load the factory program and delete any stored data, you can perform a reset. The procedure is as follows:

Turn the equipment off, press the ENTER key and hold it and start the equipment. RESET ? will appear on the screen. Confirm with the ENTER key. To cancel the action, press the ESC key. This action can be useful if you notice that the operation is not correct and before sending the equipment for technical service.

5. GUARANTEE

This product is guaranteed for one year from the date of purchase against faulty workmanship. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, S.A.U.

Any mistreatment of the equipment by unauthorised personnel not approved by J.P. SELECTA, S.A.U. cancels the guarantee automatically.