RONER

80013 (230V)
80015 (110V)

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IMMERSION THERMOSTAT

Indice

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2. GENERAL INFORMATION

- Handle the parcel with care.
- Check that contents coincide with the packing-list. If any part is damaged or missing, please advise the distributors immediately.
- Do not install or use the equipment without reading this handbook first. This instruction manual must be attached to the equipment and must be available for all users.
- If you have any doubts or enquiries, contact with ICC's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

ATTENTION!
Risk of burns

Take precautions
(use thermal gloves, clamps, etc...) to remove samples from the bath.

They could be at a very high temperature

2.1 Staff qualification

This equipment may only be used by properly trained personnel who have read and understand these instructions or have been properly trained in this equipment operation.

2.2 Use with hazardous substances

Do not treat flammable or explosive materials or those which chemically react with violence.

3. EQUIPMENT DESCRIPTION

The RONER immersion thermostats are widely used for cooking and heating baths. The unit includes a water circulation pump that allows an accurate temperature regulation and stability and improves homogeneity into the bath.

Culinary level applications:

- Cooking products previously vacuum packed (meat, fish, poultry, vegetables, terrines, pates, jams, preserves, aromatic oils, etc ...) at a controlled temperature with maximum precision and set time.
- Direct and indirect cookings are performed.
- Pasteurization (85°C) of food products cooked with traditional techniques.
- Thermal regeneration of finished food products cooked and vacuum packed.
- Ideal as a bath for heating ready food products with absolute control or those that cannot exceed a certain temperature, for example, hot siphon mayonnaise.
- When cooking under vacuum one works at a low temperature and with high precision, 58C, 62,5°C, etc ... thus achieving tenderness and juiciness.

4. PACKING LIST

The standard equipments consist on the following components:

Equipment	80013 / 80015
Instructions manual	80300

5. INSTALLATION & ASSEMBLY

Unpack the equipment and visually check that the contents are in good conditions. Keep the package for some days.

5.1 Liquid tank

A bucket or tank will be needed for the thermostat to operate. See ICC's product range on the catalogue. Our tank is made with a thermal insulation of large capacity 0.034W / mC, 40mm thick which allows for a greater stability and preservation of temperature.

Set the thermostat in the tank by a screw. Make sure that the equipment is properly secured.

Place the equipment and the tank on a flat level stable surface, suitable for the equipment weight.

5.2 Connection to the mains

Plug the equipment to a near socket suitable for the equipment power. Use the connection cable supplied or a similar one.

For your safety, the equipment must be earthed.

Before plugging the equipment, check that the voltage supply and the frequency match to the one indicated in the equipment plate.

6. STARTING UP

6.1 Control panel

1. Stop / Start cycle Key
2. Reduce value key
3. Increase value key
4. Adjustment (change between temperature / time).

6.2 Starting up

1. With the equipment installed in the bucket, fill it with water to cover the pump completely.
2. Press the start key, the display will light up and the pump will start to run.
3. Press the Up/Down keys (3/2) until the desired temperature value is selected in the green display (the above one).
4. Press key (4) if you want to select an operating time.
5. Press key (1) to start the cycle.
6. The countdown time will start counting as soon as the set temperature will be reached.
7. The heating symbol will light up and flash. 
8. Press key (1) again to stop the process. Note: during operation, temperature or time adjustment value can be changed.
9. When the time has expired, a beep will sound and the heating-element will stop heating.

REMARKS:

- * When using the time function, the green display (the above one) will indicate the remaining time. If you want to see the target temperature value, press key (4).

6.3 Safety thermostat

The equipment has a safety element to prevent overheating.

In case the thermostat triggers, press the rear button to reset the device.

6.4 Circulation pump directions


Fig 6.4.1 INTERNAL circulation.


Fig 6.4.2 EXTERNAL circulation.

The water circulation pump could operate in two configuration modes:

- INTERNAL circulation:

Water will be pumped from horizontal pump output (See fig 6.4.1).

Use this configuration to made water circulation into the bath vessel in order to homogenize the temperature from sample load.

- EXTERNAL circulation:

Water will be pumped from vertical pump output (See fig 6.4.2).

Use this configuration to made water circulation out of the bath vessel in order to temperate an external equipment.

7. MAINTENANCE & CLEANING

Electric risk

Risk of high temperature contact

MAINTENANCE:

Before operating inside, disconnect from the mains.

The manipulation of the equipment internal electronic circuits by unauthorized personnel can cause irreparable damage. Take it to one of ICC's authorized technical services.

CLEANING:

To clean the accessible parts of the equipment, use a cleaner or a similar product. Do not use aggressive solvents. To remove lime, you can use vinegar a bit lowered.

8. TECHNICAL SPECIFICATIONS

Reference code:	80013	80015
Dimensions: height/width/depth (mm)	280/180/190	
Supply voltage:	230 VAC 50 - 60 Hz	115 VAC 50 - 60Hz
Electrical power (W):	2060	1060
Weight: (Kg)	4	
External fuse:	Replaceable (Type: 10x5mm)	
Temperature range: (No refrigeration)	Ambient T ^a + 5°C ...	150°C
Display resolution:	0.1°C	

9. ERROR MESSAGES & TROUBLESHOOTING

P.FAIL Power failure. The equipment continues with the process.

Err will appear in green on the top line of the display.

(Note: if you turn off the device without pressing STOP, this error will appear when the device is restarted)

TSEG The safety thermostat is open.

SENS Open temperature probe.
The control does not detect the temperature probe.

HEAT The resistance does not work.
The equipment detects that the temperature does not increase.

Description	Cause	Resolution
It does not heat	The safety thermostat has skipped	Press the thermostat reset. See 9.1
	Heating element is melted	Contact the technical service to replace
Temperature reading does not match the external thermometer	Temperature imbalance	See 9.2 on how to act for adjustment

9.1 Reset for overtemperature

The equipment has a separate safety device to prevent that the heating element temperature rises uncontrolled, for example, when it runs out of water.

In that case, the mechanism disconnects the heating element power and, once cooled, it must be reset in order to continue using the computer.

Press the reset button on the back side of the computer.

9.2 Setting the «Offset» temperature reading

The offset (A) and the °C/°F unit (B) can be adjusted.

A. Adjust the probe reading correction (offset)

Press keys (2) and (3) continuously and at the same time (about 5 seconds). The **OFFs** message will appear.

Use keys (2) / (3) to set the actual temperature reading value. Press key (1) to confirm.

B. Presentation unit: °C / °F

Press keys (2) and (3) continuously and at the same time (about 5 seconds). The **OFFs** message will appear.

Press the key (4). The message **Unit** will appear. Using keys (2) / (3), select the unit: °C / °F. Confirm with key (1).

10. WARRANTY

This product is guaranteed for one year. It does not cover damages caused by an incorrect use or causes beyond the control of ICC.

Any manipulation of the equipment by personnel not authorized by ICC automatically cancels the guarantee.