**HORNO de MUFLA
MUFFLE FURNACE**

“SELECT-HORN-TFT”

2000376

3,6L

2000377

9L

Indice

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Index

1. General Information	12
2. Packing list	12
3. Technical specification	12
4. Information about safety in the use of the equipment	13
4.1 Risks the operator is subjected to	13
4.2 Staff qualification	13
4.3 Use with hazardous substances	13
4.4. Maintenance	13
4.5. Warranty	13
5. Equipment Description	14
6. Installation	14
6.1 Location	14
6.2 Connection to the mains	14
7. Operation	14
7.1 Inicio rápido para trabajar a una temperatura	15
7.2 SETTING screen	16
7.3 TASK screen	18
7.4 PROGRAMMING screen	19

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


1. General Information

- Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributors immediately.
- Do not install or use the equipment without reading this handbook first. This instructions manual must be attached to the equipment and must be available for all the users.
- If you have any doubts or enquiries, contact with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

2. Packing list

Product code	2000376	2000377
Size	3.6L	9L
Thermal shelf:	x1	x1

3. Technical specification

	2000376	2000377
Usable space:	3.6L	9L
Chamber dimensions: High/wide/depth (mm)	10/15/24	52/54/56
External dimensions: High/wide/depth (mm)	15/20/30	58/59/65
Power supply:	230 VAC 50Hz	
Electrical power (W):	2500	3000
Weight: (Kg)	54	70
Fuse:	Type: 10x5mm	
Temperature range:	From ... T ^a amb + 5°C ... to ... 1150°C	
Temperature resolution:	1°C	
Timed operation:	week days, hh, mm	
Temperature profiles:	Up to 10 profiles 6 segments. Segment from 1 minute to 99hour	
Data recording:	USB memory device	

4. Information about safety in the use of the equipment

Hazard risk

Electric risk

Risk of high temperature contact

4.1 Risks the operator is subjected to

- Risk of burns
- Electrical risk

4.2 Staff qualification

This equipment may only be used by properly trained personnel who have read and understood these instructions or have been properly trained in this equipment operation.

4.3 Use with hazardous substances

Do not treat flammable or explosive materials or those which chemically react with violence.

ATTENTION!

Burning risk

Take proper precautions
(use of thermal gloves, clamps,...)
to extract samples from inside the furnace.

It can be at a high temperature.

4.4. Maintenance

Unplug the equipment from the mains before removing the furnace's cover to handle inside.

The manipulation of the equipment internal electronic circuits by unauthorized personnel can cause irreparable damage. Take it to one of J.P.SELECTA's authorized technical services.

It is recommended to clean the accessible parts of the furnace with a cloth and a cleaner, or a similar solution. Do not use aggressive solvents.

4.5. Warranty

This product is guaranteed for one year. It does not cover damages caused by an incorrect use or causes beyond the control of J. P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.


Fig 1. Temperature profile

5. Equipment Description

Furnace muffles are used in sample thermal treatments, for its desiccation, cured and sterilization.

Fan assisted circulation gets a great temperature uniformity and a quick temperature rise.

6. Installation

Unpack the equipment and visually check that the contents are in good conditions. Keep the package for some days.

6.1 Location

Place the equipment on a flat level stable surface, suitable for the equipment weight.

To work in a comfortable way, please leave a free space of 50cm around the equipment and a minimum of 20cm at the rear part.

6.2 Connection to the mains

Plug the equipment to a near socket suitable for the equipment power. Use the connection cable supplied or a similar one.

For your safety, the equipment must be earthed.

Before plugging the machine, check that the voltage supply and the frequency corresponds to the one indicated in the equipment characteristics plate.

7. Operation

All the equipment functions are set with the touch screen. Gently press with your finger in the center of each icon.

When starting the computer, the starting screen will be shown (Fig 7.1)




Fig 7.1. Initial screen

- **SETTINGS:**
 - Time and date: menu to set the date and time.
 - PID and Autotuning: menu to set the parameters of temperature control.
 - Downloading recorded data.
 - Temperature correction.
 - Temperature alarm limits.
- **PROGRAMMING**
 - Allows you to program the temperature profile, by defining temperatures, times and ramps.
- **TASK**
 - Allows you to select and launch the operating temperature or a temperature profile.

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


Icons meaning:

7.2.1 Set time and date

Press on the time/date icon. Press successively on the box to select the day, month, year, minute & second adjustment.

Adjust by pressing on the icons Increase and Decrease.

7.2.2 Set parameters of temperature control and «Autotuning»

For an optimal use of the furnace, P, I, D parameters of temperature control should be the appropriate value.

With a value of the appropriate P, I, D parameters, the furnace will reach the temperature in the shortest possible time and without exceeding the set temperature. However, those values can be quite different at high and low temperatures.

The furnace comes with the following values: P=10; I=100; D=25 which are appropriate for most jobs.

These values can be modified, but it's advisable to perform an «Autotuning» when you need a better accuracy.

The «Autotuning» automatically adjusts the optimal parameters for a given temperature and a load.

Press on the PID/Autotuning icon. Press successively on the box to select the P, I, D values adjustment and modify the value by pressing on the Increase and Decrease icons.

Attention: It is not recommended to change these parameters. Use the «Autotuning» option to automatically adjust these values.

Press on the PID/Autotuning icon.

Select the «Autotuning» temperature with:  
(for example, the temperature most used)

Start «Autotuning» by pressing on the icon: 

If for some reason you want to stop, press on: 


Fig 7.2. Communication

7.2.3 Communication type

The equipment can store readings from time to time and send them to a computer or store them in a USB memory pen.

- «PORT» Set the readings storing mode:
 - RS232: Send the temperature readings to a computer by a serial cable (DB9 type).
 - USB: Store the readings in a USB memory pen drive.
- «INTERVAL» Select with the register time interval and the data sending.
- «MODE» Select the data to be stored:
 - Start: Save all temperature readings.
 - Set point: Save temperature readings, since the set point temperature has been reached.
 - PC: Send data, it is not stored.
- «SEND» Number of stored data, pending from sending.


Fig 7.2. Temperature alarms

7.2.4 Set the temperature and the temperature alarms

It is recommended to verify calibration of temperature with some regularity (for example annually) by using for instance a thermometer or a temperature pattern.

If it is necessary to correct the reading, as a result of calibration, one must enter a correction.

To do this, press on this icon:

A box indicating the current correction value appears.

Adjust with keys until you get the indicated value on the temperature pattern.

Applied correction will appear in brackets.

7.2.4.1 Temperature alarms

By pressing successively on the temperature correction box, one may enter an over-temperature and low-temperature limit. Both alarms activate a message on screen.

7.2.5 Temperature units

Press on the Temperature units icon. Press successively to change between Celsius or Farenheit degrees.

7.2.6 Set language

Press on the language icon. Press successively on the box to select among Spanish, French or English.

Adjust by pressing on the Increase and Decrease icons.

7.3 TASK screen

There are two operation modes:

- Working at a unique temperature.
- Working according to a pre-programmed temperature profile.

Press on  to change the operation mode.


Fig 7.3.1. Run screen operating at one temperature.

7.3.1 Operation at a temperature

- Select the set point temperature by pressing on the temperature indicator.
- Select a temperature increase ramp if you want to increase it in a slowly way.
The value 0°C/min makes that the ramps increases as fast as possible.
- Select an operating time. This time starts to count down when the temperature reaches the set point.
The 0°C value sets a continuous operation.
- Start by pressing on this icon: 

7.3.2 Operation with a temperature profile

- Temperature profiles must be previously programmed from the «PROGRAMMATION» display.
- Select a profile (from 0 to 9) by pressing on keys   .
- Start by pressing on the icon:  .
- In the temperature box (top-left corner) the current temperature and the set point are indicated.


Fig 7.3. Run screen operating with temperature profile.




Fig 7.3.2 Working screen.

7.4 PROGRAMMING screen

Go to PROGRAMMING screen to:

- Set the temperature profiles.
- Set the timed operation for a profile.
(For instance, you can set a time an a week day to start a profile. As : start, profile n^{er} 1 on Sunday at 19h00')

7.4.1 Programming a temperature profile

A temperature profile is a sequence of steps, each step could be a ramp or a dell.

To go from the end of a «dell» to the following «dell» could be done in two ways:

- With a «temperature ramp»
(go from 100°C to 200°C following a ramp at 2°C per minute).
- As fast as possible.

Profile parameters input:



- Select the profile number between 0 and 10 with .
- Press on the parameter to be modified. (Temp, ramp, length ...)
- Press on these icons: to enter the value.
- Modify another parameter.

Special values:

- If you enter value **OFF** in **temperature** (1st column) (you have to press Increase key) the next steps are disabled and not executed.
- If you enter **INF** in **Step length**, the furnace remains permanently at this temperature.
- If you enter **OFF** in the **ramp**, the furnace makes the temperature rise at the fastest speed possible.

Starting of temperature profile

- To set the profile starting time, press on «-:—».
- Enter time by pressing in hh:mm format.
- Press on the days of the week when the profile should be run.
- Accept the programming by pressing on .

J.P. SELECTA, s.a.u. 2015