
**TERMOSTATO BLOQUE METALICO BIO-BATH 2**  
**METALLIC BLOCK HEATER BIO-BATH 2**

**7001562**

## INDICE

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## 1 Safety

The equipment includes appropriate security measures.

Risky situations that must be respected are indicated along this manual.

### 1.1 Safety icons

The safety icons identify the risky situations and the safety measures that must be taken.

Risk of danger  
Electrical risk

Respect the indicated instructions to perform the operation.

### Electrical risks

Risk of electrical accident when accessing the areas indicated by this sign or when performing the operations indicated by this icon in this manual.

Respect the indicated instructions to perform the operation.

### Important information

- Important information for good results or for an optimum equipment operation.
- Important information to extend life of the equipment or to prevent degradation of any of the components.

### 1.2 Personnel qualification

This equipment can only be used by personnel who have been properly qualified.

This equipment can only be used by personnel who have read and understood these instructions or who has been properly qualified in the operation of this equipment.

## 2 General information

Handle the package with care. Unpack and check that the contents coincide with the "Packing List". If any part is damaged or missing, please advise the distributor immediately.

Do not install or use the equipment without previously reading this instructions manual. They are an inseparable part of the equipment and should be available to all users of the machine. Any questions can be clarified by contacting J.P. SELECTA's technical service.

Any modification, deletion or lack of maintenance of the device could cause damages. The manufacturer is not liable for the damages that may arise.

The equipment should be put in the place of low temperature, little dust, no water and no sun or strong lamp. The place should be good aeration, no corrosively gas or strong disturbing magnetic field, far away from central heating, camp stove and other hot resource. Don't place the equipment in a wet and dusty place.

### Maintenance

**Manipulation of electronic circuits by unauthorized persons may cause damages difficult to repair. Check with J.P. SELECTA's technical service.**

To clean the accessible parts, use a cloth stained with alcohol cleaner or a similar product. Do not use aggressive solvents.

## 3 Warranty

This product is guaranteed for one year. The warranty does not cover damages caused by an improper use or causes beyond J. P. SELECTA, s.a.u.


Any manipulation of the equipment by unauthorized personnel cancels the warranty benefits.


## 4 Technical specifications

**Voltage supply 100-240VAC as indicated in the machine's features plate.**

Code:

Input power	DC24V
Maximum power	40W
Average power	5 W
Temperature range	Ambient +5 - 100°C
Time range	1sec - 999sec 1min - 999min 1hour - 99hours
Accuracy of temperature	≤± 0.3°C
Display accuracy	0.1°C
Heating time (from 40°C to 100°C)	≤20min
Ambient temperature	5 - 35°C
Dimension (L×W×H)	145mm×113mm×110mm
Net weight	0.5 kg


## 5 Packing list

The standard equipment consists of the following components:

- Incubator 7001562
- Adapter
- Instructions manual 80341


## 6 Introduction

Bio-Bath 2 is a new dry bath controlled by a miniprocessor, which can be widely used for the preservation and reaction of the samples, the amplification of DNA, the pre-denaturation of the electrophoresis and blood serum coagulation.

## 7 Installation

- Place the incubator onto a level, horizontal surface.
- Insert the column connector of the adapter into the power connector of the device, and insert the other side of the adapter connector to mains power supply.
- Power on the main switch. The incubator will heat up till the setting temperature, it needs about 5min.
- Put the test tube into the block, then put on the lid.

## 8 Operation

### 8.1 Setting new procedure

Press Start/Stop continuously (more than 2 seconds) and enter into programming. Users can set procedures from p1 to p9.

----28.7----

### 8.2 Setting procedure number

Press "P" and adjust ID number through "▲" or "▼" keys.

W/  
P2 55.5 095 M  
/ \

### 8.3 Setting temperature

Press Temp, adjust temperature through "▲" or "▼" keys.

W/  
P2 55.5 095 M  
/ \

### 8.4 Setting time

Press Timer, adjust time through "▲" or "▼" keys.

When the time is "000", press "▼", then the screen changes to "KEEP", which means the instrument will run continuously. Press "▼" again and the screen displays "OFF", which means this procedure is canceled.

W/  
P2 75.5 032 M  
/ \

P2 55.5 KEEP

P2 75.5 OFF

### 8.5 Setting unit of time

Press Timer, users can set the unit of time through "▲" or "▼" when cursor is in this position.

W/  
P2 55.5 032 H  
/ \

W/  
P2 55.5 032 M  
/ \

W/  
P2 55.5 032 s  
/ \

### 8.6 Start/Stop

With the sound of "BEEP", instrument displays current temperature and heat up to 56°C.

----28.7----

P1 32.5 WAIT

"WAIT" in screen means the instrument is heating .

P1 56.0 060 M

When temperature reaches, instrument starts to count down with the sound of "BEEP"

--Finished-50.5--

When finished, press any key to enter into setting interface. Press Start/Stop continuously, it runs again. During running, press Start/Stop continuously, instrument stops running and enters into setting interface.

### 8.7 Optional setting

Press "P" and "Timer" simultaneously when power on, and the system will enter into the setting interface.

rtf	ENABLE
-----	--------

rtf	DISABLE
-----	---------

**1.RTF**

Starting RTF, the system will keep module's temperature after finishing running. Canceling RTF, the system will start to cool module's temperature after finishing running.

orf	ENABLE
-----	--------

orf	DISABLE
-----	---------

**2.ORF**

Starting ORF, if power is suddenly lost during running, the system will automatically run again the procedure after power is recovered.

Canceling ORF, the users need to select the procedure to run manually after power is recovered.

apf	ENABLE
-----	--------

apf	DISABLE
-----	---------

**3.APF**

Starting APF, the instrument will automatically run as last procedure after each boot. Canceling APF, the users need to select the procedure to run manually after each boot.

When the instrument is sold out, the APF functions are opened, RTF and ORF are closed, the users can set manually as per requirement. Both APF and ORF are opened, ORF takes precedence over APF.

Press Start/stop for a moment, the system will save the set information and exit the current interface.

Press continuously (more than 1 second), and the system won't save the set information.

### 8.8 Temperature calibration

The temperature of the instrument has been calibrated before it is sold out. But if there is deviation between the actual temperature and the displayed temperature due to some reasons, you can do as follows to correct the error.

Notes: The Instrument uses multipoint temperatures adjustment to ensure its veracity. The temperature veracity will be within  $\pm 0.3^{\circ}\text{C}$  after the multipoint temperature adjustment.

Both the circumstances and the block temperature should be lower than  $35^{\circ}\text{C}$ .

Adjustment methods as follows:

8.8.1 After the startup of the Instrument, it enters waiting interface, make sure the temperature in display is below  $35^{\circ}\text{C}$ . If the temperature is higher than  $35^{\circ}\text{C}$ , you should wait until the temperature is below  $35^{\circ}\text{C}$ .

8.8.2 Inject olefin oil into one of the cone-shaped wells, and then put a thermometer into this well (Make sure the precision of the thermometer should be within  $0.1^{\circ}\text{C}$  and the temperature ball should be absolutely immersed into the cone-shaped well). Heat insulation material is needed on the block to separate it from the circumstance. See chart A.



Chart A

8.8.3 Press Temp's "▼" and Timer's "▼" keys simultaneously, system enters calibration mode and automatically heats up to 40°C.

8.8.4 When the temperature reaches 40.0°C, the decimal place will flicker.

8.8.5 If the actual temperature of thermometer is 39.8°C. Press temp's "▲" or "▼" key, amend temperature to reach 39.8°C, then press "stop/start" key to confirm.

8.8.6 Then system will automatically heat up to next calibration. Users need to repeat above steps and input data at each calibration.

8.8.7 After 4 calibration finishes, press start/stop to confirm, screen displays "ADJ FINISHED", users can reboot instrument for normal use.

**Notes:**

During multi points temperature calibration, press "P" continuously (more than 2 seconds) to cancel the calibration.