**DRYTIME II**

**2000912**

**ESTERILIZADOR POUPINEL CALOR SECO  
POUPINEL STERILIZER**

## INFORMACIÓN GENERAL

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## GENERAL INFORMATION

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.
- 7) Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

## PACKING LIST

The standard equipment consists of the following components:

Piece	Code	Qty.
Drytime II sterilizer	2000912	1
Shelf		1
Rings		2
Heating cover shelf		1
Instructions Manual	80384	1

## TECHNICAL FEATURES

Voltage supply 115-230 V 50/60 Hz as indicated on the characteristics plate.

<b>Code</b>		<b>2000912</b>
Capacity (L)		1.7
Stability* (°C)		± 6
Internal dimensions (cm)	Height	6
	Length	24.5
	Width	11.5
Overall dimensions (cm)	Height	18
	Length	31
	Width	31
Consumption (W)		300
Weight Kg.		5

\*Measured according to the measurement procedure of J.P. SELECTA, s.a.u.

## EQUIPMENT DESCRIPTION

Poupinel DRYTIME II sterilizers are designed to allow temperature regulation from 50°C to 200°C, with a stability of ± 6°C.

It is well known that high temperatures can kill bacteria and viruses. The maximum temperature of drying time is 200°C. This device is widely used in hospitals, chemical, food, catering, and sterilizes companies, from utensils, tools to small glass, or metal items, etc.

Temperature control is done by an hydraulic thermostat, and it also includes a timer from 0 to 120 minutes for automatic disconnection.

Quick sterilization of several instruments (surgical, odontological,...) thanks to the heating by shielded elements in the base.

Built-in removable tank with extraction rings.

## INSTALLATION

Place the sterilizer on a flat horizontal, level surface, leaving a space of 15cm at the back and sides of the equipment.

Do not install the sterilizer near sunny places, heating, etc.

For an optimum operation, ambient temperature must be under 30°C.

**CAUTION!!! IMPORTANT FOR YOUR SAFETY**

**Do not use the incubator to dry or for thermal treatment of products which can give off vapours capable of making an explosive mixture.**

**Be sure that the voltage supply is the same as the one indicated on the characteristics plate.**

**Do not use the equipment if it is not earthed.**

**If you change the plug bear in mind the following :**

<b>Blue cable:</b>	<b>Neutral.</b>
<b>Brown cable:</b>	<b>Phase.</b>
<b>Yellow/green cable:</b>	<b>Earth.</b>

## OPERATION

CONTROL PANEL:

1. Main switch/timer disconnection selector.
2. Mains indicator lamp (green).
3. Timer knob.
4. Adjustment thermostat knob.
5. Heater operation indicator lamp (yellow).
6. Analog temperature reading thermometer.

## OPERATION

1. Switch the equipment on to the correct power supply. Please check the product power at the back before usage.
2. Place the material to be disinfected inside the machine.
3. Adjust the time. The red pilot will light off when disinfection is complete. Usually, sterilizing the goods at 125°C for 15 minutes could be enough to kill 99.9% bacteria and viruses.
4. Adjust the temperature with the left hand knob. It could be from 0 to 200°C. The green pilot will light on meanwhile the oven is operating.
5. Make sure the door is close before turning the timer button.
6. Keep children away from the machine.
7. Switch it out when you don't use it to avoid possible fire.

NOTE: When the sterilizer is tempered, it is possible to make the equipment work directly in a temporized way for consecutive sterilization, because in these conditions the rise temperature time is almost nil. When the user works with the same load at the same temperature, the time of temperature rise is known. This permits the starting up of timer from the first moment; nevertheless, the sum of rising and sterilization time must be prefixed.

### COMMENTS:

- \* The graduation of the thermostat scale is only a guide as the real temperature inside the oven is the one indicated in a calibrated external thermometer.
- \* It is not advisable to work the sterilizer at temperatures below 100°C as the thermostat does not regulate correctly in these conditions.
- \* All regulation and readjustment of temperature must always be made on an upward scale, that is before choosing any nominal value, the control must be turned completely to the left of the scale. A more exact regulation of temperature is obtained in this way as it compensates the inevitable movement of the thermostat spindle.

## **MAINTENANCE**

**Before removing the casing, disconnect the equipment from the mains.**

**The manipulation of the internal electronic circuits of the incubator by unauthorized personnel can cause irreparable damage. Take it to one of the J.P.SELECTA, s.a.u. authorized technical services.**

CLEANING:

For the cleaning of the different parts use the following products:

Cleaning of stainless steel: alcohol.

Cleaning of plastic: alcohol with cotton duster.

## **SPARE PARTS**

**To guarantee the safety of the equipment, all spare parts must be purchased from J.P.SELECTA, s.a.u.**

## **GUARANTEE**

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.