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## 1. General Information

- Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- Do not install or use the equipment without reading this handbook first. This instruction manual must be attached to the equipment and must be available for all users.
- If you have any doubts or enquiries, contact with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

## 2. Packing List

Model	3000982
Shaker with support trays	1
Trays	9
100 graphics pack (ref.0231740)	1
Instructions manual	80390

## 3. Technical Specifications

Bag capacity (litres)	54
Useful measures: height/width/depth (cm):	385/445/390
Displacement length (mm):	40
Electric power (W):	310W
Supply voltage	110-230VAC 50/60Hz
External fuse	Replaceable (type: 10x5mm)
Temperature range	5°C ... 60°C
Display resolution	0.1°C
Weight (Kg):	26
Sound level (dB):	40
Loading and unloading programmable time:	99s
Loading and unloading break counter:	99
Fail counter by supply failure:	99
Timed operation	week, hour, minute
Date and time indication:	hh/mm/ss dd/mm/yy
Temperature profiles	Up to 10 profiles of 6 segments. Segment from 1m to 99 hours
Parameters protection with password:	Fixed password 1111 (shaker) 0000 (incubator)
Thermo log, Temp. margin	From -10°C to 40°C      1 lap = 1 week

## 4. Equipment Description

This equipment consists of two equipments that can work independently.

### 4.1 Shaker

The stirrer consists of a fixed base on which rests a movable platform that allows a lateral swinging movement. A fixed structure that supports the sliding trays is installed on the platform.

### 4.2 Incubator

The chamber is made from methacrylate that allows to view inside the chamber during incubation. Its fan circulation system ensures the homogeneity of temperature. The temperature regulation is made by an electronic circuit board.

The incubator is used in thermal treatment of samples for drying, curing, cultivation or preservation at a constant temperature.

Heating and cooling system by semiconductors and temperature regulation is electronical with digital reading. The cooling system is affected by ambient temperature. The minimum temperature increases with the ambient temperature according to the following table:

Ambient temperature	Minimum temperature
25°C	5°C
26°C	6°C
27°C	7°C
28°C	8°C
29°C	9°C
30°C	10°C

The incubator allows programming temperature profiles to optimize the heat treatment. It has an operation log system for later downloading data to a computer. In addition, it is equipped with a paper thermograph that works independently of the incubator. It has a USB front connection to download the logs to memory.

### 4.3 Paper chart recorder

The incubator includes a mechanical thermal recorder that reads the temperature inside the chamber. It works with an AA type battery so that if you do not want the temperature to be recorded, it is necessary to remove the battery.

## 5. Safety in the use of the equipment

Hazard  
risks

Electrical  
risks

Risk of high  
temperature contact

### 5.1 Risks the operator is subjected to


- Risk of burns
- Electrical risks

### 5.2 Staff qualification

This equipment may only be used by properly trained personnel who have read and understood these instructions or have been properly trained in this equipment operation.

### 5.3 Use with hazardous substances

Do not treat with flammable or explosive materials or with those which chemically react with violence.

## 6. Installation and Start-up

### 6.1 Packaging

Once the equipment is unpacked, visually check that it is in good conditions. Keep the package for some time just in case you will need it.

### 6.2 Location

Place the equipment on a flat, level, stable and adequate surface to the equipment's weight, keeping a space of about 25cm from the wall, to allow ventilation and movement of the equipment.

Level the equipment if necessary, by using the four adjustable level rubber bumpers on the stirrer base.

Please bear in mind the equipment weight already loaded, to determine the weight the table must support.

### 6.3 Connections to the mains

Before plugging the equipment, check that the voltage supply and the frequency corresponds to the one indicated in the equipment characteristics plate.

For your safety, the equipment must be earthed.

If you change the plug, bear in mind the following:

- |                     |          |
|---------------------|----------|
| Blue cable:         | Neutral. |
| Brown cable:        | Phase.   |
| Yellow/Green cable: | Earth.   |

The audible alarm of power failure is activated during the first connection of the shaker to the mains.

For a correct operation, the equipment connection voltage must be  $\pm 10\%$  of nominal voltage. To avoid problems and failures, it is important that the equipment works within these margins.

The incubator open door detector that works on the platelet agitator must be connected to it by means of the cable from the back side, as shown in the pictures.

Under these conditions, when the incubator door is opened, agitation is interrupted to facilitate loading and unloading work.

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Fig. 6.4 MAIN SCREEN (START)

### 6.4 Starting

Choose an electrical outlet near the equipment and appropriate for the equipment power. Use the supplied connection cable or a similar one.

At the back of the incubator there is an additional power outlet to connect the agitator power and thus be able to turn off the assembly only through the incubator main switch.

It starts with the **MAIN SCREEN** when pressing the main switch (Fig. 6.4). **RUN & STOP** buttons are activated according to the state they were when switching the incubator off with the O/I switch the last time.

From this screen, you can access the different operating options:

- Settings menu. (See chapter 8.5)
- Simple work mode / programming. (See chapter 7)
- Start / stop the incubator.

10:42:20 24/01/18

The status bar indicates in which menu is located, the date and time and if there is any type of communication connected (in the example, the USB connection icon appears).

### 6.5 Installation of the paper chart recorder

1. Open the recorder lid by using the key
2. Remove the needle
3. Unscrew the shaft handle and remove the paper
4. Pull out the shaft to remove the recorder
5. Insert the battery into the rear part
6. Insert the recorder again
7. Put the needle again and close the recorder
8. Close the recorder lid again

Use a good quality alkaline battery to ensure a duration of more than 6 months.

Note: the recorder needle is supplied with a cap that protects it during transport. Remove it before using the recorder.

## 7. Shaker Operation

The platelet stirrer has been designed for a quiet and continuous blood bag operation.

To enable a continuous operation monitoring, the equipment records all break times, both due to loading and unloading, as unforeseen ones, due to power failure.

### 7.1 Starting-up

Start the equipment by pressing the main switch.

By a (short) press of the button, it stops for loading and unloading.

The stop lasts for a determined «T» time. After this time, the stirrer starts automatically.

With a long press you can access the record menu.

### 7.2 Break counter

One can access the break counter by the single button.

This counter can be reset to zero.

The break lasts for a determined «T» time. After this time, the stirrer starts automatically.

### 7.3 Power supply fails

The computer registers the power failures.

For every stop, the computer registers date, time and duration of the stop.

Registration of stops can be reset.

### 7.4 Password

Access the above parameters is password protected.

Password is 1111.

Password can be momentarily unlocked, to facilitate access to counters and registers before initiating treatment.


### 7.5 Settings


The agitation speed is factory set. If it is necessary to modify this speed, remove the blue rear plug and adjust the speed by using a flat screwdriver.

### 7.6 Parameters selection rotatory button

This equipment has an only button in the control panel, which has a couple of movements: circular and rotatory.

## 7.7 Menus Diagram

MAIN SCREEN (START)


Fig 8.1 Screen to set temperature and time.

Fig 8.2 Screen to stop a cycle.

## 8. Incubator Operation

All the equipment functions are established through the touch screen. Press gently with your finger in the center of each icon.

This screen shows the status of the oven before disconnecting it by means of the main switch (0/I).

This allows the oven to resume its task in the event of a power cut.

### 8.1 Set temperature and working time

For most applications, this is the usual work method. If you need to create temperature profiles or program the start-up of the oven automatically, follow the steps in section 7.3.

To select the working temperature:

- From the **MAIN SCREEN**, click on the value to be modified. The frame of the modifiable parameter will be boxed in yellow. In the example, the temperature can be modified.

The **+** and **-** keys appear to increase / decrease the value of the adjustment.

- After a few seconds without pressing any key, the yellow frame disappears and the adjusted value is selected.

- The same operation can be performed to modify the working time or the temperature ramp.

A time of **--h --m** establishes a continuous operation.


### 8.2 Start / Stop the cycle

- From the **MAIN SCREEN** press:

**RUN** to start the cycle (the icon will change its appearance with the **STOP** message).

**STOP** to stop the cycle (the icon will change its appearance with the **RUN** message).

Note: When the duty cycle is in operation, a heat or cold symbol appears depending on the type of energy that the heating element is generating.

MAIN SCREEN (START)


Fig. 8.3 Temperature profiles

Fig. 8.4 Editing temperature profiles

Fig. 8.4.1 Modify temperature, ramp and time parameters.

Fig. 8.4.2 Modify maximum and minimum temperature limits.

## 8.3 Programming temperature profiles

A temperature profile is a sequence of steps defined by a temperature and a time. It also allows you to define a schedule for the start / end of the process.

The equipment allows defining up to 10 profiles of 6 steps each.

From the main screen, press this button to work in "Programming" mode; the temperature profile will appear graphically as shown in Fig 8.3.

To start the program, simply press the **RUN** key (the button will change its appearance with the **STOP** message).

If during the program operation you wish to stop the process, press the **STOP** button.

From screen 8.3, click on this button to return to the main screen.

## 8.4 Editing temperature profiles

From screen 8.3, press the **EDIT** button to select and modify the program. The screen will show the appearance of Fig. 8.4.

Press the **EDIT** button to modify the temperature, time and ramp parameters, as shown in Fig 8.4.1.

Press the **ALARMS** button to define the ranges between which the system should work. If the reading of the temperature is above or below the defined values, a warning message will appear telling that the programmed limits have been exceeded. These alarms are only informative, they do not stop the operation of the equipment.

Press the **CYCLICAL** button to repeat indefinitely the selected program.

### 8.4.1 Set temperature, time and ramp parameters

The **+** and **-** keys allow you to add steps to the profile and the left / right arrows allow you to move between the different phases of the profile.

Press the **BACK** button to return to the main editing screen of profile 8.4.

### 8.4.2 Set maximum and minimum temperature limit values

If you want the oven to show an alarm message when the temperature is outside margins, set a value for over temperature and for low temperature.

Press on the value and use the Up/Down keys to adjust these limits.

Press the **BACK** key to return to the main editing screen of profile 8.4.

Fig 8.4.3 Selection date and time of program start

### 8.4.3 Selection of date and time of program start

From screen 8.4.3 you can plan the week days and the start time of the cycle.

Press "Activate" to activate programming.

Press the **BACK** button to return to the main editing screen of profile 8.4.

**NOTE:** If a program is running, the cycle will not turn off. It has to be stopped.

Fig. 8.4.4 Selection of working program

### 8.4.4 Selection and saving working program

To select one of the 10 programs, press the **SEL** button from screen Fig. 8.4 to select the working program of the oven.

Press the **SAVE** button to save all the parameter changes established in the program.

MAIN SCREEN (START)

### 8.5 Settings

Press this button from the main screen to enter the settings screen as shown in Fig 8.5

Modifying the settings parameters can affect the equipment operation. Read carefully the entire chapter so that the equipment works according to the process needs.

Most sensitive temperature control parameters are protected by a password. See the corresponding chapter.

Parameters settings are grouped according to their function as shown below:

#### LANGUAGE

Allows you to select the screen language.

#### DATE/TIME

Adjusts date and time of the oven. This adjustment is important if you work with programming.

#### UNITS

Sets parameters units.

#### BRIGHTNESS

Sets brightness of the screen backlight.

#### ENTER PASSWORD

Enters the password for accessing the calibration menu.

#### CHANGE PASSWORD

Changes initial factory password.

#### ALARMS

Sets margins out of which the oven gives an error.

#### LOGS DOWNLOAD

Press this button and the control will allow downloading of data stored in memory to a USB memory connected in the front part.

#### INFORMATION

Shows the firmware version and gives access to the menu for its update.

Fig. 8.5.1 Language selection

### 8.5.1 Language

When entering this setting, screen Fig 8.5.1 appears. Use the **Up/Down** arrows to select the language you want to appear in all menus.

Press the **BACK** key to return to the settings menu or the **HOME** key to return to the initial screen.

Fig. 8.5.2 Date / Time selection

### 8.5.2 Date/Time

Screen Fig 8.5.2 appears when you enter this setting. Click on the parameter to be modified and use the **Up/Down** arrows to select the corresponding value.

Press the **BACK** key to return to the settings menu or the **HOME** key to return to the initial screen.


Fig. 8.5.3 Units selection

### 8.5.3 Units

When entering this setting, screen Fig 8.5.3 appears. Use the **Up/Down** arrows to select the parameter from which you want to change the unit.

Fig. 8.5.3.1 Temperature unit selection

Fig 8.5.3.1 shows the screen for selecting the temperature units. By using the **Up/Down** arrows you select between °C and °F. If this value is changed, all displayed values will be referenced to this unit.

Press the **BACK** key to return to the settings menu or the **HOME** key to return to the initial screen.

Fig. 8.5.4 Setting the screen brightness

### 8.5.4 Brightness

When entering this setting, screen Fig 8.5.4 appears. Use the **Up/Down** arrows to select the brightness percentage of the screen.

Press the **BACK** key to return to the settings menu or the **HOME** key to return to the initial screen.


Fig. 8.5.5 Setting alarm margins


Fig. 8.5.6 Password entry

Fig. 8.5.6.1 Temperature Offser modification


Fig. 8.5.7 Password modification

Fig. 8.5.8 Logs download

## 8.5.5 Alarms

When entering this setting, screen Fig 8.5.5 appears. Use the **Up/Down** arrows to select the value of over temperature or low temperature you want for the temperature range.

This setting is simply a warning, it does not stop the operating cycle. This notice will be reflected in the memory log.

Press the **BACK** key to return to the settings menu or the **HOME** key to return to the initial screen.

## 8.5.6 Temperature Offset Adjustment

The Offset adjustment is protected by a password to prevent any user from modifying this parameter.

 Press the **CAL** button to access this setting.

If the **CAL** button is pressed without entering the password, it will indicate that a valid password value is necessary to continue with the action. Therefore, it is necessary to select **ENTER PASSWORD** option before pressing the **CAL** button.

The default password value for the calibrator is “**0000**”.

Press on the temperature Offset parameter and use the **Up/Down** keys to modify this value.

Press the **BACK** key to return to the settings menu or the **HOME** key to return to the initial screen.

The password value remains valid until the oven is switched off and on again. Perform this action to prevent another user from modifying the value.

## 8.5.7 Password Change

If you want to change the factory password, use this menu and follow the instructions that appear to change the password.

Press the **BACK** key to return to the settings menu or the **HOME** key to return to the initial screen.

## 8.5.8 Logs Download

To download the data stored in the control, select the option **LOGS DOWNLOAD**. Screen Fig. 8.5.8 will appear.

Connect a memory pen drive to the front connector (4) and follow the instructions on the screen.

There are two download options: either download all logs or the last 10 only.

During the loading of data in the USB memory, information about % of loading appears.

Press the **BACK** key to return to the settings menu or the **HOME** key to return to the initial screen.


Fig. 8.5.9 Information screen

Fig. 8.5.9.1 Firmware update

### 8.5.9 Information. Firmware update.

The equipment is tested and checked before shipment but there may be improvements in the control operation.

If it is necessary to update the firmware version, press the **INFORMATION** option and Fig 8.5.9 appears.

Press the **UPGR** key and follow the instructions that appear on the screen. A screen similar to Fig 8.5.9.1 will appear.

To obtain the file with the update, it is necessary to know the **DEVICE ID** code that appears in screen Fig 8.5.9. Load this file into a **USB** memory and insert it into the front connector.

Attention: during the update process DO NOT turn off the equipment or remove the **USB** memory from the connector. The equipment could run out of control.

The update process lasts a few minutes. During the update process the control is turned off several times, but this is normal.

Press the **BACK** key to return to the settings menu or the **HOME** key to return to the initial screen.

## 9. Maintenance

Before cleaning, disconnect the computer from the mains.

Do not spill liquids on the equipment.

### CLEANING:

External parts of the equipment can be cleaned with a little alcohol and water solution. To clean the different parts, it is advised to use:

Alcohol for cleaning of stainless steel.

Alcohol with cotton or a soft piece of cloth for cleaning plastic.

## 10. Warranty

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. Selecta, s.a.u

Any manipulation of the equipment by personnel not authorized by J.P. Selecta, s.a.u. automatically cancels the warranty.