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# “T4000”

## 4120824

# Contenido

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## **INTRODUCTION**

- Careful handling to avoid collision and falling instruments in use.
- Avoid contact with water or other liquids during use.
- Do not put the instrument in the sunshine for a long time. After use, it should be stored in a cool, dry and ventilated place.
- If you do not use the instrument for a long time, you should unplug the power supply to avoid accidents.
- This instrument is not suitable for use in harsh environment, high temperature, low temperature or strong magnetic field interference, which may lead to instrument damage.
- If there is any problem with the instrument, please contact the dealer or the company. Do not disassemble the instrument by yourself. If disassembled, the company will no longer be responsible for the warranty.

## **GENERAL FEATURES**

- Colour LCD display
- Intelligent menu operation
- Multiple automatic calibration
- Differential signal measurement mode, stable and reliable
- Manual and automatic temperature compensation
- Two relay control switches
- High & low alarm and hysteresis control
- 4-20mA & RS485, Multiple output modes
- Multi parameter display simultaneously shows – pH/ORP, Temp, current, etc.
- Password protection to prevent misoperation by non-staff.

## **SAFETY IN THE USE OF THE EQUIPMENT**

Before installing or operating this equipment, you should read this manual and be aware of potential hazards to avoid them.

Upon receiving the equipment, carefully open the package, check whether the equipment and accessories are damaged by transportation and whether the accessories are complete. If you find any damages, please report it to JP SELECTA immediately and keep the package for return.

This is an analytical measuring and control instrument with a high level precision. Only trained or authorized persons should carry out the installation, configuration and operation of the equipment. Make sure the power cable is physically disconnected from the power supply when connecting or repairing.

When handling the equipment, make sure it is turned off and unplugged. Although it has been designed for safety operations, any equipment can involve potentially dangerous processes or procedures. The main potential risks associated with this particular equipment are:

Misuse injuries

Electric shock injuries

Potential danger is located more in the measurement points where the electrodes are installed than in the equipment itself, given the applications where they are installed may be in contact with potentially dangerous substances.

All precautions must be maintained and potential hazards must be considered. It is dangerous to misuse this equipment or to ignore the instructions, rules and warnings.

Do not exceed the maximum specified operating conditions, disconnect the equipment from the electrical network before handling or accessing the terminal strip.

## TECHNICAL SPECIFICATIONS

Measurement range:	pH: -2 ~ 16 pH; ORP: -1999- ~ + 1999 mV.
Unit:	PH, mV.
Resolution:	PH: 0,01pH; ORP: 1 mV
Reading error:	PH: +/- 0.02pH; ORP: +/- 2mV
Temperature:	0 ~ 150.0°C (depends on the sensor)
Temperature resolution:	0.1°C
Temperature accuracy:	+/- 0.3°C
Temperature compensation:	0 ~ 150.0°C (manual or automatic)
Stability:	pH: <0.02pH / 24h; ORP: <2mV / 24h.
Current outputs:	Two: 4 ~ 20mA, 20 ~ 4mA, 0 ~ 20MA.
Communication Output:	RS485 MODBUS RTU.
Relay control outputs:	Two, 3A 250VAC, 3A 30VDC.
Power supply:	85 ~ 265VAC. Optional: 9 ~ 36VDC, consumption 3W
Working conditions:	no magnetic field interference except geomagnetic field
Working temperature:	-10 ~ 60°C
Relative humidity:	<90%
Protection degree:	IP65
Weight:	0.6kg
Dimensions:	98x98x130 mm.
Installation opening size:	93x93 mm
Installation modes:	Panel and wall mounting

## INSTALLATION

The equipment can be installed inserted inside a cabinet, by making a 93\*93mm hole. Or on the wall, by using the stainless steel support, which comes with the kit, and fixing it to the wall with M4 plugs and screws.

### Connectors:

N/C- No connection	NO1- High setting relay
N/C- No connection	COM1 Setting relay High common
Temperature sensor	NC1 Idle high setting relay
Temperature sensor	NO2 Low setting relay
Ent. Signal +	COM 2 Adjustment relay Common bass
Ent. Signal - (Ground)	NC2 Low setting relay standstill
N/C	
N/C	Phase- (85 ~ 265V) Live wire
N/C	AC N Neutral wire
Cur 1+: S. Current Channel 1	PE Earth cable

COM - Common Current Output

Cur 2+: S. Current Channel 2+

RS485 B Communication

RS485 A Communication

## **Electric connection**

Be careful when unpacking to avoid damage to the equipment. Make sure that all services are compatible with the equipment.

Remove the back cover from the equipment by loosening the screws. Through the wall bushing, insert the signal cables.

Connect the power supply, input and output signals, relay contacts and equipment motherboard. The sensor cable length is usually 5-10 meters. Cable ends are inserted into the connectors and tightened securely. Always use reliable, cutting-edge connections, do not improvise.

Place the back cover, tighten the four screws firmly to maintain IP65 degree of protection. Tighten the cable glands of the signal cables.

Make sure the equipment is grounded and connected to an electrical outlet with the correct voltage.

The power supply must incorporate an Earth-leakage circuit breaker (ELCB) or a Residual current circuit breaker (RCCB) to protect the user from electric shock in case of misuse or accident.

Potential hazards should always be the first consideration when deciding on a suitable location for equipment. Leave enough space around the equipment.

## STARTING UP

Although our equipment has been designed to be installed and configured very easily, we recommend to be carried out and reviewed by competent personnel. Connect the equipment to the power supply. This screen will appear:


On the main screen, in big digits, you can see the measurement variable reading of pH or Orp, in the picture above, it would be pH.

On the bottom left there is the pH measurement in mV. And on the right hand side, the temperature measurement.


There are two LEDs under the screen on the left hand side: SP1 that indicates the status of the high alarm relay (Blue off. Red on) and SP2 that indicates the status of the low alarm relay. (Blue off. Red on). In addition, there are four blue round buttons (from left to right):

- Menu / ESC: to enter the programming menu and return to the previous one. If pressed down, it returns to the measurement screen.
- Right Arrow: used to move through the options menu.
- Down Arrow: to move through the options menu and input values.
- ENT: to confirm the options. If we press the (M/E) button, the following programming menu screen will appear:


With the key (Down Arrow) we place ourselves on the different menus and pressing (ENT) we open the selected menu, for example: "Configure".

With the key (Down Arrow), we select between Temperature Sensor. By pressing (ENT), we open the sensor menu.

With the (Down Arrow) key, select Type or Unit. Press (ENT) to open the Sensor menu.

PH or ORP appears with the (Down Arrow) key. We select one and validate it with the (ENT) key.

Press the (M / E) key and return to the immediately previous menu. Thus we will configure all the operating parameters of the equipment.

If we keep the (M / E) key pressed, we will return to the main screen.

## Programming menu structure

The following is the menu structure of this equipment.

Press the (Menu) key to enter menu setting mode:

- If you are asked to enter the calibration security password, press the (Down) or (Right) key to set it, and then press the (ENT) key to confirm the calibration security password.

There is no initial password, enter directly by pressing the (ENT) key.

## CONFIGURATION

<b>Sensor</b>	Type: PH, ORP Unit: pH, mV
<b>Temperature</b>	S. Temp: NIC2. kQ, NIC10 kQ, Pt100, Pt1000 Temp offset: 0,0000 Temp input: automatic, manual Temp. unit: °C, °F

### CALIBRATION:

<b>Standard Cal.</b>	USA: 7.00, 4.01, 10.01: Id. auto. pH 7.00 first NIST: 6.86, 4.01, 9.18: Id. auto. pH 6.86 first ORP Sol. standard:   235mV (default) Offset 1, Slope 1, Offset 2, Slope 2 correction
<b>Field Cal.</b>	Cal. Field, Aj. Offset1, Aj. Slope1, Aj. Offset2, Aj. Slope2

### ALARM:

<b>Relay 1</b>	Estado: On, off High / low set point: high alarm, low alarm Limit value Hysteresis
<b>Relay 2</b>	Estado: On, off High / low set point: high alarm, low alarm Limit value Hysteresis

### OUTPUT:

<b>Current 1</b>	- Channel: Main, temperature - Output option: 4-20mA, 0-20mA, 20-4mA - Upper limit - Lower limit
------------------	---

Current 2	- Channel: Main, temperature - Output option: 4-20mA, 0-20mA, 20-4mA - Upper limit - Lower limit
RS 485	- Baud rate: 4800, 9600, 19200BPS
Parity:	None, odd, even
Stop Bit:	1 Bit, 2 Bit
	Network node

SYSTEM:

Language	- chinese, english
Display	- Display speed: low, standard, medium, high - Backlight: Saving, Bright
Software	- Software version: 19-1.0 - Password: 0000 (default) - Serial num: SHCY 2002 00466
Default values	- 1. No, 2. Yes
Programming Current Outputs	-1: 4-20, 20-4 -2: 4-20,20-4
Relay test	- Relay 1. Relay 2

## CALIBRATION

Press the (MENU) key to enter setup mode and select the "**Calibration**" menu. The following menu will appear:

CALIBRATION: <b>Standard</b>	USA: 7.00-4,01-10,01. (Automatic)
	NIST: 6,86-4,01-9,18 (Automatic)
	ORP Standard: 235mV (Default)
<b>Field</b>	Field calibration, Offset, Slope.

## Standard solution calibration

Select Standard Solution Calibration. There are two options, you can choose between:

EE. UU: 7.00, 4.01, 10.01

NIST: 6.86, 4.01, 9.18,

Once the selection is made, press the (ENT) key to confirm and enter the standard solution calibration mode.

If the equipment has been calibrated, the display shows the status of the calibration. Press the (ENT) key again to recalibrate the equipment if needed.

If you are asked to enter the calibration security password on the screen, enter the password value using the (Down) and (Right) keys. The default will be 0000. Then press the (ENT) key to confirm the calibration security password.

The following display will appear::

PH calibration: after entering calibration mode, the equipment is displayed as above.

The equipment automatically identifies the standard liquid. It first calibrates the midpoint (example 7, .00pH), then calibrates 4.01 pH or 10.01 pH. The corresponding mV value will be displayed on the left side of the screen.

When calibration is complete, the offset and slope will be displayed on the right side of the screen.

Press the (Menu) key after marking the calibration points, and exit directly.

During the calibration process, an error message will appear on the screen when the standard fluid is wrong.

Calibration results: the slope of the glass electrode (> 0.90) is allowed and that of the metal antimony electrode (> 0.80) is allowed.

## **ORP calibration**

Press the (Menu) key to enter the setup menu. Select the ORP standard liquid calibration and enter the known standard liquid value (default 235m), press the (ENT) key to enter the standard solution calibration mode.

If the equipment has been calibrated, the display shows the status of the calibration. Press the (ENT) key if you need a re-calibration.

After entering the calibration mode, the corresponding mV value will be displayed on the left side of the screen.

When calibration is complete, the slope will be displayed on the right side of the screen. The pH value of the buffer solution was measured at 25 ° C.

To calibrate the instrument using an automatic identification buffer, you will need a standard pH buffer that matches any of these values. Before using automatic calibration, select the correct buffer table (see "Buffer Table").

Before calibration, the sensor must be kept immersed in the Calibration Standard to ensure the stability and accuracy of the calibration and reading values.

## **Field calibration**

Select the on-site calibration methods: Linear Calibration, Offset Adjustment, or Linear Adjustment.

### Linear Calibration

When the known pH or ORP value is entered in this field, the equipment will automatically correct the data.

### Offset Adjustment

Compare the data on the lab or laptop with the data on the computer. If there are errors, the data can be modified by this function.

### Linear Adjustment

Linear values after "field calibration" are stored in this item, factory default is 1.00.

## **COMPONENTS INCLUDED**

Our equipments are supplied complete and ready for basic operation.

## **PACKING LIST**

The equipment must be carefully unpacked and the components must be verified with the packing list.

Any omission or breakage must be reported to JP. SELECTA s.a.u.

<b><u>DESCRIPTION</u></b>	<b><u>QUANT</u></b>
T4000 equipment	1
Bulwarks Pg10	2
Side fixing clips	2
Wall mounting bracket	1
User manual	1

## **MAINTENANCE**

Daily maintenance:

In general, this equipment does not need daily maintenance. If there is a fault, you can contact our company and make adjustments and repairs with the help of our technical service.

Simple verification method:

Set the zero potential parameter (Set) 0.0 mV and the slope (Slope) to 1.00, and then short-circuit the meter terminals "PH +" and "PH-", and the meter should show 7.00 or be within range of basic mistake.

Connect a 2K resistor to the "29" TEMP + and "30" / TEMP terminals and the meter temperature should show around 27.5°C or be within the basic error range.

If the display is not correct, it means the meter is not working. You will need to contact the manufacturer for maintenance.

For the maintenance of the electrode, attention should be paid to the following points:

- (1) When the electrode is immersed in a recirculation circuit, the speed of the incoming water sample must be kept as constant as possible.
- (2) The electrode head is provided with a protective cap containing liquid from the electrode to maintain the moisture of the membrane. Remove the protective cap before use and put it on when not in use.
- (3) A new electrode should be immersed in a 3.3 MKCL solution for 24 hours to ensure the accuracy of the calibration.
- (4) Calibration is recommended prior to installation. When installing the measurement point, the glass bulb of the composite electrode must not collide to avoid damage.
- (5) When the measurement point can be dry for a long time, the protective cap should be placed on the electrode, filled with a concentration solution of 3.3 mkc2. It should not be dried, otherwise the electrode could be damaged.
- (6) If there is dirt in the glass chamber on the top of the electrode, it can be cleaned with 0.1n hydrochloric acid, then activated by dipping it in a 3.3mkcl solution.

**Do not touch the electrode by hand.**

To ensure correct measurement for a long time, the electrode must be cleaned, maintained and calibrated regularly.

The antimony electrode will produce oxidation after being used for a period of time, and the electrode head will turn black, fine sandpaper can be used to polish the electrode head and remove the oxidation part, to ensure the sensitivity of electrode measurement.

(7) Keep the electrode cable and connector clean and dry, not wet or with water.

(8) The electrode will be replaced if it fails or is damaged. It is recommended to replace it every year.

When the electrode is not in use, it should be soaked in the maintenance solution.

For the preparation of the immersion solution:

Take a pack of pH4 buffer, dissolve it in 250ml pure water, and then add 56 g of pure analytical KCl, heat it properly and stir until completely dissolved.

(note: the electrode is easy to be spoiled when stored in dry state)