



Flake ice makers

Refrigerated cabinets “Medilow” XS, S, M, L & LG

Refrigerated medical cabinets, S, L & M

Upright freezers “Templow” S, M, L, P and Z

Storage refrigerators “Stocklow” LS & GS

Chest freezers “Combatemp”

page 71

page 72

page 73

page 74

page 75

page 75



LOW TEMPERATURE

Model selection table

QUICK OVERVIEW

Equipment type	Model	Forced Air	Regulation	CAPACITY			
Refrigerated Cabinets				105 litres			
+2 ... +8 °C	Medilow XS	Yes	Analogical	2101286			
Cabinets				180 litres	236 litres	544 litres	544 litres glass
+1 ... +40 °C	Medilow-S-M-L-LG	Yes	Digital	2101270	2101271	2101272	2101273
Storage refrigerators				236 litres	352 litres	330 litres	
+0 ... +15 °C	Pharmalow-S-L	Yes	Analogical	2101274	2101275	-	
+2 ... +8°C	Pharmalow-M	Yes	Digital	-	-	2101287	
Storage freezers				513 litres			
-14 ... -25 °C	Tempow-M	No	Digital	2101279			
				257 litres			
-15 ... -32 °C	Tempow-S	No	Digital	2101278			
				600 litres			
-10 ... -35 °C	Tempow-L	Yes	Digital	2101280			
				289 litres (Two compartments)			
+2 ... -40 °C	Tempow-YC	Yes	Digital	2101288			
Storage refrigerators				670 litres	670 litres Glass		
+1 ... +8 °C	Stocklow LS	Yes	Digital	2101283	-		
+1 ... +8 °C	Stocklow GS	Yes	Digital	-	2101284		
-17 ... -22 °C	Stocklow LSC	Yes	Digital	2101289			
Chest Freezers				206 litres	215 litres	343 litres	352 litres
-14 ... -28 °C	Combatep	No	Digital	2101040	-	2101041	-
-10 ... -40 °C	Combatep	No	Digital	-	2101181	-	2101182
Flake ice maker				Production 58 Kg	Production 89 Kg	Production 142 Kg	
	"COM"	Compressor	-	5020104	5020105	5020106	



Flake ice maker "COM-50", "COM-85" and "COM-150"

"COM-50" MODEL: PRODUCTION 45 kg-24 HOURS.
"COM-85" MODEL: PRODUCTION 85 kg-24 HOURS.
"COM-150" MODEL: PRODUCTION 150 kg-24 HOURS.



MODELS	Part No.	24h.(kg) Production	Tank Capacity (Kg)	Height / Width /Depth cm			Power W	Weight Kg
EF-50	5020104	58	20	79	47	60	470	58
EF-85	5020105	89	20	79	47	60	550	58
EF-150	5020106	142	60	136	51	55	900	75

Features: See page. 190.





Refrigerated cabinet “Medilow-XS”

FOR TEMPERATURES FROM 2°C TO 8°C.
IDEAL FOR SMALL SPACES.

FEATURES

Epoxy covered external case.
Interior ABS molded, resistant to chemicals with LED light built-in.
Internal fan for temperature homogenization.
Supports for adjustable in height shelves.

CONTROL PANEL

Analog temperature controller.
On/Off switch internal light.
On/Off switch external light.

MODELO

Código	Capacity litres	Height / Width / Depth (interior) cm			Height / Width / Depth (exterior)cm			Power W	Power HP	Nº Shelves	Weight Kg
2101286	105	146,5	28	27	182,5	36	38,5	110	1/8	8	65

Supplied with 8 steel trays, PVC laminated. For additional tray. Part No. 1001861



Refrigerated cabinets “Medilow” S, L & LG

WITH FORCED AIR FAN CIRCULATION.
WITH DIGITAL ELECTRONIC TEMPERATURE CONTROL.
TEMPERATURE RANGE CONTROLLABLE FROM +2 °C TO +40 °C. WITH A STABILITY OF ±1.5 °C. RESOLUTION 1 °C.
CAPACITY: 180, 610 AND 544 LITRES.

SAFETY: CONFORMS TO EN61010.1. OVER TEMPERATURE SAFETY THERMOSTAT INSTALLED.



Model L.



Model Medilow LG, double glass door.



Model Medilow S. Capacity: DBO 2 blocks

APPLICATIONS

Determination of BOD at 20°C, Measurement of enzyme activity at 25°C. Storage and Preservation of samples, Incubation, Cosmetics, Botany, Food Products, etc.

TEMPERATURE CONTROL

Electronic temperature control with digital display.
Pt 100 temperature probe.
Probe failure alarm.

FEATURES

Exterior cabinet covered in epoxy and interior enclosure-molded in ABS, resistant and easy to clean.
Hermetic and vibration-free compressor group.
Air circulation turbine inside the enclosure.
Side hole for entry of connections, probes, etc.
Power socket outlets inside (power maximum 500W).
Level adjustable front foot supports.
Change of door opening direction.
Lockable door.
Door can be hung to open left or right.
Door lock.

CONTROL PANEL

Temperature control.
Live-time temperature display.
Pre-selected temperature display.
Increase value, push button.
Decrease value, push button.
Program parameter selection, push button.
Configure push button.

MODELS	Part No.	Capacity litres	Door glass	Height / Width / Depht (interior) cm			Height / Width / Depht (exterior) cm			Number of shelves	Internal power points	Power W	Weight Kg
MEDILOW-S	2101270	180	-	70	51	44	85	60	85	3	2	250	36
MEDILOW-L	2101301	610	-	145	60	56	164	75	95	3	2	330	73
MEDILOW-LG*	2101302	544	SI	145	60	56	164	75	95	4	2	330	82

Comes complete with steel shelves, PVC laminated.

*Anti mist double glazed insulated door, with LED interior light.

SPARE PARTS Shelves

For Medilow S/L/LG Part No.	2101270	2101301	2101302
Shelves	1001813	1001986	1001987



Refrigerated cabinet “Pharmalow” L

WITH FAN ASSISTED CIRCULATION.
TEMPERATURES RANGE FROM +1 °C UP TO +15 °C.



APPLICATIONS

Products that require a constant temperature below ambient, preservation of pharmaceuticals, medicines, microtubes, pharmaceuticals, medicines, microtubes, etc.

Automatic defrosting with water evaporation.
Level adjustable foot supports.
Door direction change.
Door with lock.

COMMON FEATURES

Exterior cabinet and door in light gray lacquered metal and interior enclosure in white sheet metal with non-toxic finish.
Metallic door and interior led light.
Hermetic compressor unit.
Air circulation turbine throughout the interior enclosure for greater homogeneity.

CONTROL PANEL

Temperature regulator thermostat inside.
Digital thermometer temperature reader.
Digital electronic temperature controller.

SPARE PARTS Shelves

For Pharmalow /L

Shelves **1001985**



MODEL

Part N	Capacity liters	Height / Width / Depht (interior) cm	Height / Width / Depht exterior (cm)	Power W	Number of shelves	Weight KG
PharmalowL 2101300	327	145 47 44	169 60 65	180	4	58

Comes complete with steel shelves, PVC laminated.

ACCESSORIES Need to be factory fitted.

OPTIONAL

Drawer with 3 compartments and 12 partitions.

Code **1001315**



Storage refrigerator with double glazed anti mist door “Pharmalow” M

DIGITAL ELECTRONIC TEMPERATURE CONTROL FROM +2 TO +8

MODEL

Part No.	Capacity litres	Height / Width / Depht (interior) cm	Alto / Ancho / Fondo (exterior)cm	Power W	Power HP	Number of shelves	Weight Kg
2101287	330	145 54 47	184,5 62 60	575	1/4	6	113

Comes complete with 6 steel shelves, PVC laminated.

ACCESSORY “Pharmalow” S - L and M. Must be factory fitted.

Circular graphic register independently powered.

Temperature range: -10 °C up to 40 °C. Precision $\pm 2\%$ at 12 °C.

Chart speed selection of 24 hours or 7 days. Disc Ø 125mm

Graphic register dimensions:

18 cm high x 22 cm wide x 35 cm deep.

Powered by an LR-6 1.5V battery. Battery life 1 year.

Part No. **2101621**



Model Pharmalow M with 6 shelves.

SPARE PARTS

Part No.

1001860 Set 2 additional guides + shelf



Upright freezers “Templow” S, M,L and YC



APPLICATIONS

Preservation and storage of products that require to be preserved at low temperatures.

Upright freezer “Templow” S

DIGITAL TEMPERATURE CONTROL,

FIXED TEMPERATURE: -15 °C, -18 °C, -21 °C, -25 °C AND -32 °C.

8 DRAWERS WITH INDIVIDUAL FREEZER COMPARTMENTS.

FEATURES

Epoxy covered metallic external case.
Strong, molded, chemically resistant ABS interior.
Manual defrost with water evaporation system.
Door can be hung to open left or right.

CONTROL PANEL

Mains power ON/ OFF switch.
Temperature thermostat regulator.
Temperature indicator lamp.
Audible and visual alarms.
Push button alarm reset/stop.
Push button with quick freeze indicator.



Upright freezer “Templow” M

TEMPERATURE RANGE FROM -14 °C TO -25 °C.

7 INDIVIDUAL FREEZER COMPARTMENTS WITH 14 BASKETS.

FEATURES

Epoxy covered external case.
Strong, chemically resistant, ABS molded case.
Specially coated “Rilsan”, wire mesh baskets.
Manual defrost, with water evaporation system.

The door can be hung to open left or right.
Lockable door.

CONTROL PANEL

Temperature thermostat regulator.
Digital temperature display.



Upright freezer “Templow” L

WITH FAN ASSISTED AIR CIRCULATION.

TEMPERATURE RANGE FROM -10 °C UP TO -35 °C.

COMES WITH 4 SHELVES.

FEATURES

External and internal case made completely of
AISI 304 stainless steel, with water sump and drain.
Chemically resistant.
Acclimatized hermetically sealed CFC free compressor.
Automatic defrost with electrical water evaporation
system.
Regulable shelves.

Adjustable front level supports.
The door can be hung to open, left or right.
Door lock.

CONTROL PANEL

Mains power ON/ OFF switch.
Electronic temperature regulator.
Digital temperature display.



Refrigerated cabinet with double compressor “Templow” YC

UPPER CASE: TEMPERATURES FROM 2°C UP TO 8°C. CAPACITY: 189 LITRES.

LOWER CASE: TEMPERATURES FROM -10°C UP TO -40°C. CAPACITY: 100 LITRES.

3 SHELVES AND 1 DRAWER IN UPPER CASE AND 3 DRAWERS IN LOWER CASE.

CARACTERÍSTICAS

External case covered in epoxy.
Stainless steel internal easy for cleaning, with also
stainless steel wire basket.
Double compressor that offers maximum stability and
constancy of temperature.
Gradual adjustment of the trays for a versatile use of
the interior.
Door with lock and key.
USB output for data recording.

CONTROL PANEL.

Digital electronic control.
Alarms of low and high temperature, door open and
power failure.
Backup battery for alarms.
Upper case temperature indicator display.
Lower case temperature indicator display.
Parameters value increase/decrease.



NEW

MODELS	Part No.	Capacity litres	Range °C	Height / Width / Depht (interior) cm	Height / Width / Depht (exterior) cm	Power W	Weight Kg	
TEMPLOW S	2101278	257	-15 to -32	151 44 42	184 60 63	110	81	Comes with 8 drawers
TEMPLOW M	2101279	513	-14 to -25	151 61 56	186 75 75	345	109	*Comes with 14 baskets
TEMPLOW L	2101280	600	-10 to -35	155 51 65	215 70 83	600	130	*Comes with 4 shelve
TEMPLOW YC	2101288	289	**2 to -40	- - -	183 70 64	360	136	Comes with 3 shelves and 4 drawers

*Comes with steel PVC laminated shelves and baskets. **in different compartments



Storage refrigerators “Stocklow” LS, GS and LSC

WITH FAN ASSISTED CIRCULATION.

MODEL LS, WITH METAL DOOR. TEMPERATURE CONTROLLABLE FROM +1 °C TO +8 °C.

MODEL GS, WITH DOUBLE GLAZED ANTI MIST DOOR. TEMPERATURE CONTROLLABLE FROM +1 °C TO +8 °C.

MODEL LSC, WITH METAL DOOR. TEMPERATURE CONTROLLABLE FROM -17 °C TO -22 °C.



NEW
DESIGN

APPLICATIONS

Preservation and storage of products that require below ambient temperatures.

FEATURES

Chemically resistant.

Tropical hermetic compressor.

Automatic defrost with electrical water evaporation system.

Support rack and shelves.

Adjustable front level supports (Models LS and GS).

Supplied with wheels (Model LSC).

The door can be hung to open left or right.

Door lock.

There are two versions:

Stocklow LS and LSC which is completely made of metal, both interior and exterior, AISI 304 stainless steel door.

Stocklow GS inside made of AISI 304 stainless steel and epoxy white external case with glazed door that enables internal view.

CONTROL PANEL

Mains power, ON/OFF switch.

Electronic temperature regulator.

Digital temperature display.

Interior light only (Double glazed door model).



Model Stocklow LS and LSC, with metallic door.



Model Stocklow GS, with metallic door.

SPARE PARTS

Código **1001858** Set 2 guides + additional tray

MODELS	Part No.	Capacity litres	Door	Height / Width / Depht (interior) cm			Height / Width / Depht (exterior) cm			Number of shelves	Power W	Weight Kg
STOCKLOW LS	2101283	600	Metallic	154	61	63	213	73	84	3	250	120
STOCKLOW GS	2101284	600	Glass	154	61	63	213	73	84	3	250	125
STOCKLOW LSC	2101289	600	Metallic	154	61	63	213	73	84	3	500	125

Comes complete with steel shelves, PVC laminated.



Chest Freezers “Combatep”

CONTROLLABLE TEMPERATURE RANGE FROM -14 °C TO -28 °C.

CONTROLLABLE TEMPERATURE RANGE FROM -10 °C TO -40 °C.

APPLICATIONS

Suitable for preservation and freezing of sample.

FEATURES

Epoxy covered external case. Interior of aluminium, with guide supports for baskets or trays.

The interior of the door is made of easy to clean, strong, chemically resistant, moulded ABS.

Drop down lockable door with compensating support hinges. Hermetically sealed compressor, silent and vibration free.

Automatic interior light when the door is opened. (all models except Combatep S-30.)

Digital electronic temperature control. N.T.C. Sensor.

CONTROL PANEL

Temperature control.

Digital live and preset temperature display.

Push-button temperature selector.



MODELS	Part No.	Capacity litres	Freezing capacity in 24hrs Kg	Temperature up to °C	Height / Width / Depht (interior) cm			Height / Width / Depht (exterior) cm			Power W	Weight Kg
COMBATEMP S-30	2101040	206	28	-28	74	63	53	92	75	71	160	38
COMBATEMP L-30	2101041	343	38	-28	74	100	53	92	113	71	200	56
COMBATEMP S-40	2101181	215	30	-40	63	89	41	92	113	76	230	57
COMBATEMP L-40	2101182	352	40	-40	65	117	50	92	137	81	400	72

Note: Part No. 2101040 is supplied with 1 basket. Part No. 2101041 is supplied with 2 baskets. Part No. 2101181 is supplied with 3 baskets and Part No. 2101182 is supplied with 2 baskets.

SPARE PARTS Baskets

For Combatep. Part No.	2101040 / 2101041	2101181	2101182
Baskets	1001810	1001811	1001812